

User Inspired Designs Since 1990

Push and Transfer Gloves

Open finger design enables independent donning and doffing by persons with limited hand function


2 3/4" of dual Velcro closure for a secure fit

Leather or Non-Slip Rubber Surface enhances grip for transfers, push and slow downs. Protects palms and wrists against skin abrasions

Open fingers, half thumb, wrap-around design

Exterior made of 100% Leather

Padded in palms and wrist



ADI's Push and Transfer Gloves are one size fits most. The Velcro closure allows for size adjustment Available in black All Leather or black Leather with Rubber Grip