

AirLogic™ Anterior Trunk Supports*: Slim-Cut

AEL EXCLUSIVE! AirLogic – AEL exclusive patented technology – is a series of breathable fabrics and structural layers that maximize air movement under and through the durable posture support. Air Mesh enables perspiration to pass through the mesh keeping users cool and dry. The durable Stretch panels incorporate a structured flexibility that prevents tearing.



Slim-Cut Stretch



Slim-Cut Non-Stretch



Slim-Cut Zippered

AirLogic™ Slim-Cut Anterior Trunk Support

Stretch & Non-Stretch



Stretch Part No	Non-Stretch Part No	Size	A	B	C	D	E	F	Webbing Width
14256	14286	X Small	8 3/4"	14 1/12"	5"	8 1/4"	1 1/4"	2"	1"
14257	14287	Small	9 1/4"	17"	6"	10"	1 1/4"	2"	1"
14258	14288	Medium	10 1/4"	20"	7 1/2"	11 1/4"	1 1/4"	2 1/4"	1"
14259	14289	Large	11 1/4"	22"	8 1/4"	12 1/2"	1 1/2"	2 1/2"	1"
14260	14290	X Large	12 1/4"	24"	9"	13 1/2"	1 1/2"	2 3/4"	1"

AirLogic™ Slim-Cut Anterior Trunk Support

Zippered

Stretch Part No	Size	A	B	C	D	E	F	Webbing Width
14158	Medium	10 1/2"	21"	7"	11 3/4"	2 1/2"	2 1/4"	1"
14159	Large	11 1/2"	23"	7 3/4"	13"	2 3/4"	2 1/2"	1"
14160	X Large	12 1/2"	25"	8 3/4"	14"	3"	2 3/4"	1"

To size Anterior Trunk Supports, match dimensions in the table with the illustration



Suggested For Physical Conditions:

- Trunk or Shoulder Protraction
- Muscle Weakness/Atrophy
- Spasticity
- Scoliosis
- Trunk Lordosis
- Trunk Kyphosis Causing Forward Trunk Flexion
- Trunk or Shoulder Rotation
- Body Temperature Control

*AirLogic is an exclusive patented technology of AEL