

Belts | Hip Belts | Easy Ordering

Hip Belts – Easy Ordering

Hip Belt with Push-Button Buckle



Part No	Push-Button Buckle Style	Webbing Width	Width per Pad (2 Pads/Belt)	Pull Style
14633	Non-Spring Metal	1"	6"	Single-Pull
14643	Spring	1 1/2"	8"	Single-Pull
14644	Spring	1 1/2"	10"	Single-Pull
14663	Non-Spring	1 1/2"	8"	Single-Pull
14653	Spring	2"	10"	Single-Pull
14654	Spring	2"	12"	Single-Pull
14645	Spring	1 1/2"	8"	Dual-Pull



Hip Belt with Side-Release Buckle

Part No	Webbing Width	Width per Pad (2 Pads/Belt)
14603	1"	6"
14604	1"	8"
14623	2"	10"

Suggested For Physical Conditions:
 Posterior Pelvic Tilt • Pelvic Obliquity • Low Tone

Suggested For Wheelchair Conditions:
 Positions the Pelvis in the Chair

Hip Stabilizing Belts – Easy Ordering


4-Point Hip Belt with Push-Button (Spring) Buckle



Part No	Size	Webbing Width	Width per Pad (2 Pads/Belt)
14749	Small	1"	5 1/2"
14750	Medium	1 1/2"	7 1/2"
14751	Large	2"	9 1/4"

For optimum effectiveness, center the buckle and pad between the hips. Use the top straps to pull the hips back, stabilizing the lower spine. The bottom straps should be pulled down across the pelvis to position it at 90°. Once straps are installed and pad is centered, use the pull rings to adjust the length.

4-Point Hip Belt with Side-Release Buckle



Part No	Size	Webbing Width	Width per Pad (2 Pads/Belt)
14754	Small	1"	5 1/2"
14755	Medium	1 1/2"	7 1/2"
14756	Large	2"	9 1/4"

Suggested For Physical Conditions:
 Anterior Pelvic Tilt • Pelvic Rotation • High Tone

Suggested For Wheelchair Conditions:
 Positions the Pelvis in the Chair

Hip Belts – Itemized Ordering

Hip Belts and Hip Stabilizing Belts are now available itemized with separate part numbers for the Hip Belt and Hip Pads. If you prefer to itemize your Hip Belt and Pads, just follow these simple steps:

1. Using the tables below, choose the style of Hip Belt or Hip Stabilizing Belt
2. Then, select the size
3. Use BOTH part numbers listed in the table: one part number for the Hip Belt plus a part number for the Hip Belt Pad Set

Note: Hip Belt and Hip Belt Pad part numbers must be ordered together

Hip Belt with Push-Button Buckle

ORDER BELT

&

ORDER PAD

Hip Belt					Hip Belt Pad (Set)
Part No.	Push-Button Buckle Style	Webbing Width	Pull Style	Pad Size	Part No.
14633-99	Non-Spring Metal	1"	Single-Pull	6"	9999-01
14643-99	Spring	1 1/2"	Single-Pull	8"	9999-01
14644-99	Spring	1 1/2"	Single-Pull	10"	9999-01
14663-99	Non-Spring	1 1/2"	Single-Pull	8"	9999-01
14653-99	Spring	2"	Single-Pull	10"	9999-01
14654-99	Spring	2"	Single-Pull	12"	9999-01
14645-99	Spring	1 1/2"	Dual-Pull	8"	9999-01

Hip Belt with Side-Release Buckle

Part No.	Webbing Width	Pull Style	Pad Size	Part No.
14603-99	1"	Single-Pull	6"	9999-01
14604-99	1"	Single-Pull	8"	9999-01
14623-99	2"	Single-Pull	10"	9999-01

Hip Stabilizing Belts – Itemized Ordering

Hip Stabilizing Belt with Push-Button Buckle

ORDER BELT

&

ORDER PAD

Hip Belt				Hip Belt Pad (Set)
Part No.	Size	Webbing Width	Pad Size	Part No.
14749-99	Small	1"	5 1/2"	9999-02
14750-99	Medium	1 1/2"	7 1/2"	9999-02
14751-99	Large	2"	9 1/4"	9999-02

Hip Stabilizing Belt with Side-Release Buckle

Part No.	Size	Webbing Width	Pad Size	Part No.
14754-99	Small	1"	5 1/2"	9999-02
14755-99	Medium	1 1/2"	7 1/2"	9999-02
14756-99	Large	2"	9 1/4"	9999-02

Belt Options

Strap Guides


Keeps the posture support straps from slipping off the shoulder or from rubbing the neck of the client

Part No	Webbing Width	Unit of Measure
14581	3/4" & 1"	1 set (2 ea)
14582	1 1/2"	1 set (2 ea)



Anchor Band (Frame Mounting Hardware)

Part No	Tube Diameter	Unit of Measure
80187	7/8"	1 set (2 ea)
80188	1"	1 set (2 ea)



Replacement Strap Kit

Includes four straps with male buckles at one end and belt anchors at the other

Part No	Belt Style	Webbing Size	Unit of Measure
14840	Dynaform®	3/4"	1 set
14841	Dynaform®	1"	1 set
14842	AirLogic™	1"	1 set

Belt Modifications

Non-Stretch Modification

Convert the panel of a Stretch Dynaform to a Non-Stretch panel. The Non-Stretch style is a durable padded vinyl that is easy to wipe clean

Part No	Compatible With
14000	Any Hip Belt, Hip Stabilizing Belt, Dynaform Chest Support, Dynaform Y-Style, Dynaform Upper-Body, or Dynaform Backpack Style Support

Shoulder Positioners

POSAInc® Shoulder Positioners

- Swings up and out of the way by pivoting to the sides
- Locking knob prevents accidental lifting
- Multi-positioning design puts pads in perfect spot
- Metal-base is custom-formable
- Pads are covered in Neoprene



Part No	Pad Shape	Pad Width	Pad Length	Unit of Measure*
18405	Small	1 1/2"	6"	1 ea
18406	Medium	2"	8"	1 ea

*Specify the client's side if ordering an uneven quantity

Suggested For Physical Conditions:

Designed for rigid support to help position the client's shoulders and upper body. The multi-positioning posts and joints allow for the pad to be positioned for optimum support.