

# e·motion® e·fix® Compatibility Chart

The following chart represents a sample of compatible wheelchairs.  
For a wheelchair not shown, please contact us for more information.

Difficulty level of installation is indicated to the right of each wheelchair.  
E = Easy M=Moderate D=Difficult

## Invacare

Tracer	EX, EX2, SX, IV	Yes	E
	Insignia	Yes	E
	MyOn	Yes	E
9000	XT#, XDT#, SL#	Yes#	E
9000	XT Recliner	NO	
Action	MVP, MVP Junior	Yes	E
Action	Xtra	Yes	E
Action	Patriot	Yes	E
Action	Pro Spin X4	Yes	E
Action	A4* / A4 Ti*	Yes*	E
Action	Solara	Yes	E
Action	Solara 2G/3G (E.fix Battery Tray Required)	Yes	M
Top End	Reveal °Δ≠	Yes°Δ≠	D
Top End	Terminator/SS*Δ (With Zero degree fixed camber)	Yes*Δ	M
Top End	Terminator/SS*Δ (With reversible camber)	Yes*Δ	E
Top End	Crossfire+*Δ / TiT6 +*Δ / T6A+*Δ	Yes+*Δ	M
Top End	T7 +*Δ (Note: Chair must be built with T6 Camber System)	Yes+*Δ	M

## Sunrise Medical

Quickie	2, HP, HP-TI	Yes	E
Quickie	2 Lite *	Yes*	M
Quickie	GPV, GPV Competition	Yes	E
Quickie	TS, Zippie TS ≠	Yes≠	D
Quickie	GPS, GP, GP Titanium	Yes	E
Breezy	All Series	Yes	E
Quickie	Kidz	Yes	E
Quickie	Zippy, Zippy 2	Yes	E
Quickie	Zippy GS	Yes	E
Quickie	Zippy Zone	NO	E
Quickie	Revolution	Yes	E
Quickie	Titanium (TI) +*	NO+*	
Quickie	XTR	Yes	M
Quickie	LX	Yes	E
Quickie	LXI	Yes	E
Quickie	Iris ≠ (E.fix Battery Tray Required)	Yes≠	E
Quickie	M6 ∞	Yes∞	E
Quickie	GTX	Yes	M
Quickie	GT/GTi + ∞ ∞	Yes+	D
Quickie	Q7	NO	
Quickie	QXI Δ≠	YesΔ≠	D

## TiLite+

AERO R	Yes	M
EVO, ZRA, AERO Z	Yes	M
YR, YRA, YG, YGS	Yes	M
SX, X,	Yes	E
TR, TRA, TRC	Yes	M
2GX, AERO X	Yes	E
ZR	Yes	M
Z10, TX, XC	NO	

## Freedom Designs

	Freedom 2,	Yes	D∞
Gizmo	Tilt-in-Space ✓	Yes	M ✓
	SP3	Yes	M

## Ki Mobility

Catalyst 4	Yes	E
Catalyst 5, 5Ti, 5VX, Spark	Yes	E
Tsunami, Little Wave, Rogue **≠	Yes**≠	M

## Colours

Eclipse	Yes	E
Zephyr, Razorblade,	Yes	D
Spazz, Spazz G	Yes	D
Saber	Yes	D

## LifeStand ≠

LS	Yes	M
LSE	Yes	M
LSA	Yes	M

## Pride Mobility

LiteStream XF	Yes	E
---------------	-----	---

## Levo ✓≠

Levo	All Chairs	Yes	M
------	------------	-----	---

Due to sudden manufacturer product changes, this list is subject to change without notice.

One-Arm drive systems are not compatible with e.motion.

\*FMS Anti-Tip Supports are NOT compatible. Must use wheelchair brand Anti-Tip Supports. °Compatible only with FMS anti-tippers.

+Zero Degree Camber ONLY. ✓Additional support for frame required. ∞Restrictions may apply. #Space Saver arm rests only.

∞E.fix only installation must be done at FMS. ≠Installations must be done at FMS. ΔCompatible only with E.motion