

Dream Series

Name	Dealer	Reference Name
Phone	Ship Address	
email	City, State, Zip	

Standard side models


Our standard side (single) door system comes 20" above the mattress and has doors on both sides.

Twin 36x76 Full 54x74


- 1000 2000 Standard side door system, fixed (non-adjustable) bunkie board
- 1500 2500 Standard side door system, fixed height, manual adjustable head
- 1500-D Standard side door system, fixed height, manual head and foot
- 1800-S Standard Side door system, Semi electric (80" long)
- 1800 Standard side door system, Hi Lo, Full Articulation (electric)
- Extended Twin 80" long




Select Door Panels




Clear view windows



Clear view windows with air flow panel



Solid panel



All airflow panel

High side models

Our high side (bifold) door system comes 32" above the mattress with door on one side standard.


Twin 36x76 Full 54x74

- 1200 2200 High side (bi-fold) door system, fixed (non-adjustable) bunkie board
- 1700 2700 High side ((bi-fold) door system, fixed height, manual adjustable head
- 1700-D High side (bi-fold) door system, fixed height, manual head and foot
- 1900-S High side (bi-fold) door system, semi electric (80" long)
- 1900 High side (bi-fold) door system, Hi Lo, Full Articulation (electric)
- Extended Twin 80" long



Select Door Panels


Top door




Clear view windows



Clear view windows with air flow panel




Solid panel



All airflow panel


Bottom door




Clear view windows



Clear view windows with air flow panel



Solid panel



All airflow panel

- Dual side open (doors on both sides) for high side Part #900

Dream Series Quote /Order Form

Wood* Maple Oak

Stain* White paint (upcharge*) Natural Fruitwood
 Michael's Cherry Asbury Rosewood (upcharge*)

Head/Footboard shape:*

Colonial



Hemisphere



Straight
(required for canopy)



Accents*

Chantilly



Soothing Wave



Serenity



Custom engraving:
(write exactly what you want engraved)

Corner Posts:

Squared off on top



Turned finial
(not compatible with canopy)

Part #203



Feet: Leveling nickel glides 3" Casters (wheels) Part #120

Transfer height:

24 inch (standard) 21.5 inch (Standard side only) Part #510 30 inch Part #500

Do you need clearance for a floor lift? NO YES 5 1/2" YES 6 1/2"

Dream Series Quote /Order Form

Accessories:

NOTE: it is often extremely difficult or impossible to add or remove many of these accessories after the bed is built because many are structural changes.

- | | |
|--|-------------------|
| <input type="checkbox"/> IV Pole and mounting bracket | Part #150 |
| <input type="checkbox"/> Storage Drawers under bed rail | Part #450 |
| NOTE:
Drawers must have 30" transfer height and CANNOT be used with electric deck | |
| <input type="checkbox"/> Variable height option for fixed models | Part #455 |
| <input type="checkbox"/> Clear view window in headboard | Part #825 / 827 |
| <input type="checkbox"/> Clear view window in footboard | Part #826 / 828 |
| <input type="checkbox"/> Airflow panel in headboard | Part #825a / 827a |
| <input type="checkbox"/> Airflow panel in footboard | Part #826a / 828a |
| <input type="checkbox"/> Clear Back Panel for Hi Side beds | Part #800 |
| <input type="checkbox"/> Airflow Back Panel for Hi Side beds | Part #800a |
| <input type="checkbox"/> Foam pad lining inside of the bed (Twin) | Part #700 |
| <input type="checkbox"/> Foam pad lining inside of the bed (Double) | Part #710 |




NOTE: once pads are added the user must ALWAYS have pads or unsafe gap will occur. Similarly, pads cannot be added without extensive frame modification after bed is built because frame is oversized to accommodate. Plan accordingly.


COLOR for foam pad cover (see color chart) _____

- | | |
|--|-----------|
| <input type="checkbox"/> Access Port (for feeding tubes) | Part #725 |
|--|-----------|

Draw circles on outlines below for desired port location(s)



Headboard
Inside View



Footboard
Inside View

Dream Series Quote /Order Form

Mesh Safety Enclosure (Canopy) Option:

- Mesh Canopy Top (Twin) Part#255
- Mesh Canopy Top (Double) Part#256



Canopy Color _____

Heavy duty mesh Yes

For pad cover and mesh safety enclosure (canopy)

Fabric Colors



Polyester - Urethane Coated
Washable, mildew and water repellent

All images on this page are for general reference only. Electronic reproductions should not be used for exact color match. In addition, slight color variations will occur with any dyed fabric. We reserve the right to make product changes without notice.

availability of color
is not guaranteed

Updated January
2019



Comments / Special Requests
