



Safety Beds for Special Needs Kids

Haven Dimensions

36 x 76 Mattress (Std. Twin) - No Casters

36 x 80 Mattress (Ext. Twin) - No Casters

	3000-LO	3000	3500	3500-D
Transfer Height <i>(floor to top of mattress)</i>	13.75	20.75	20.75	20.75
Clearance Under Side Rail	0	4	4	4
Inside Height <i>(top of mattress to top of canopy)</i>	49.5	49.5	49.5	49.5
Outside Height <i>(floor to top of wooden exterior frame)</i>	63.5	70.5	70.5	70.5
Outside Depth <i>(front opening to back)</i>	42.75	42.75	42.75	42.75
Outside Length <i>(left to right)</i>	82.75	82.75	82.75	82.75

	3000-LO	3000	3500	3500-D
	13.75	20.75	20.75	20.75
	0	4	4	4
	49.5	49.5	49.5	49.5
	63.5	70.5	70.5	70.5
	42.75	42.75	42.75	42.75
	86.75	86.75	86.75	86.75

Casters - Add 3" to Transfer & Outside Heights

***ALL DIMENSIONS IN INCHES**