

# BESTMOVE™

## STA400 | STANDING TRANSFER AID

The BestMove™ STA400 is the perfect blend of form and function for patient transfers. With the new TS30760 Raising belt, users and caregivers will be able to easily transfer the patient from a sitting position to a standing as well as from a standing position to a sitting position.




- **400lb** capacity
- Shin brace for secure transports
- First choice for any fall prevention program
- Support strap available



# | STA400



CLOSED SWIVEL SEATS




OPEN SWIVEL SEATS

## | DIMENSIONS


### Key dimensions of STA400

maximum capacity	400Lb / 182Kg
overall height	43.2 in
overall base length	29.9 in
overall width	24.5
minimum seat height	26 in
maximum seat height	33.5
front-end base height	3.6
external front-end leg width	16.5
maximum kneepad height	21
horizontal handlebar width	20
net weight	60lb/27kg
shipping weight	73lb/33kg


\* SPECIFICATIONS PROVIDED ARE A GENERAL GUIDELINE  
AND ARE SUBJECT TO CHANGE WITHOUT NOTICE



FRONT VIEW OF STA400



SIDE VIEW OF STA400



TOP VIEW OF STA400