


PUSH-BUTTON

SIZE	PART NUMBER
S38	N/A
M40	HB409-M40
M46	HB409-M46


Sub-ASIS COMPATIBLE BELT

M46	LB409-M46
-----	-----------



PLASTIC SIDE-RELEASE

SIZE	PART NUMBER
S38	HB419-S38
M40	HB419-M40
M46	HB419-M46




REHAB LATCH

SIZE	PART NUMBER
S38	N/A
M40	HB429-M40
M46	HB429-M46

Sub-ASIS COMPATIBLE BELT

M46	LB429-M46
-----	-----------



X-SMALL METAL PUSH-BUTTON

SIZE	PART NUMBER
S38	HB439-S38
M40	N/A
M46	N/A

-A1 + **FLAT-MOUNT™**
END FITTING PART NUMBERS
S (1") HKA1-R-1
M (1-1/2") HKA1-G-1
L (2") HKA1-B-1

-B1 +

-A2 + **CINCH-MOUNT™**
END FITTING PART NUMBERS
S (1") HKA2-R-1
M (1-1/2") HKA2-G-1
L (2") HKA2-B-1

-B2

note

description
HIP BELT, 4-POINT, REAR-PULL, TECHNICAL SPECS

drawing no. D16-1116-32	date originated 11/16/2016	drawn by D. DOUGAL	rev. B
-----------------------------------	-------------------------------	-----------------------	------------------