

NEW!

Bodypoint Essentials H-Style Shoulder Harness

Our updated H-Style Shoulder Harness provides strong but cost-effective shoulder positioning for improved head and trunk control. This new, value-engineered design offers a cleaner appearance and a more ergonomic shape for greater comfort and performance, along with the industry-leading quality you have come to love and expect from Bodypoint.

CORDURA® fabric outer layer and proprietary edge binding eliminate wrinkles, preserve skin integrity and comfort, and increase durability and dirt-resistance

Non-stretch, closed-cell foam for superior resiliency and pressure distribution, plus soft, wrinkle-free fabric lining for optimum comfort

D-ring pivot point allows flexibility in lower strap angle without compromising fit; New pad contours improve comfort over lower ribs and minimize puckering when strap angle changes

Clean appearance with no top-stitched webbing

Adjustable sternum strap with a side-release buckle and easy-grab folded end stop allow for a fine-tuned fit



Rear-pull top straps have folded end stops with teal accents for easy adjustment and side-release buckles for quick removal

Part Numbers

Size	Part Number
Small	SH224S-B5
Medium	SH224M-B5
Large	SH224L-B5
X-Large	SH224XL-B5

Replaces the discontinued part numbers SH210 and SH220.

End-Fittings

Includes four high-strength Plastic Flat-Mount™ end-fittings, now reinforced with embedded stainless steel washers. Can also be mounted with Cam Buckles or Cinch-Mount™ end-fittings (sold separately).



Plastic Flat-Mount™

- B5

ISO Standards For a Strong Position

The Essentials H-Style Harness is tested to the same rigorous standards as all Bodypoint products, guaranteeing performance and long-term durability backed by our limited lifetime warranty.

For more about our quality management system, certified as conforming to ISO9001:2015, see bodypoint.com/standards.

CAUTION: Whenever an anterior trunk support is used, a properly-adjusted pelvic support should also be worn to prevent sliding.

