





**FRONT SIDE (AEROMESH)**



9mm AEROMESH FABRIC FACE PAD

**BACK SIDE**



VELCRO HOOK PATCH

VELCRO-COMPATIBLE LOOP FABRIC

PART NUMBER	SIZE
SP103S	SMALL
SP103M	MEDIUM
SP103L	LARGE

note

description

**AEROMESH CALF PANEL, TECHNICAL SPECS**


drawing no.  
D11-0831-33

date originated  
8/31/2011

drawn by  
D. HAGAN

rev.

**A**



note

description

**AEROMESH CALF PANEL, TECHNICAL SPECS**

drawing no.  
**D11-0831-33**

date originated  
8/31/2011

drawn by  
D. HAGAN

rev.  
**A**