

TOTAL LENGTH

S38	127cm (50")
M46	152cm (60")
L62	165cm (65")

ADJUSTMENT STRAP RANGE

S38	7.5cm (3")
M46	10cm (4")
L62	11.5cm (4-1/2")

BUCKLE

CLOSED CELL FOAM PAD, 7mm THICK

WEBBING WIDTH

S38	2.5cm (1")
M46	4cm (1-1/2")
L62	5cm (2")


DOUBLE-FOLD  
END STOP

PADDED LENGTH

S38	38cm (15")
M46	46cm (18")
L62	62cm (24-1/2")


PAD WIDTH

S38	5cm (2")
M46	5.5cm (2-1/4")
L62	7cm (2-3/4")



PLASTIC  
SIDE-RELEASE

SIZE	PART NUMBER
S38	GB215-S38-B5
M46	GB215-M46-B5
L62	GB215-L62-B5



DURALATCH


SIZE	PART NUMBER
M46	GB285-M46-B5
L62	GB285-L62-B5

INCLUDES PLASTIC FLAT-MOUNT™ END FITTINGS

S (1") HKA31-25MM-1

M (1-1/2") HKA31-38MM-1

L (2") HKA31-50MM-1



description

ESSENTIALS HIP BELT, 2-POINT, CENTER-PULL, TECHNICAL SPECS

drawing no.

D21-0721-01

date originated

7/21/2021

drawn by

M. KOSH

rev.

A