

NEW!

Grooved Mushroom Joystick Handle


A revolutionary new joystick handle, the Grooved Mushroom takes the best features of several previous Bodypoint® models and advances them to the next level. Its hand-filling shape, top and side grooves, and height-adapability make it ultra versatile to meet a wide range of user needs.



Part Numbers

PC110A	For 4.8 mm (3/16") stem size
PC110B	For 6.4 mm (1/4") stem size

Easily cut to desired height with a utility knife



60 mm
(2-3/8")

28 mm
(1-1/8")