

# Ensure all-day joystick access

## 2 SELECT MOUNTING CLAMP

### Tri-Lock™ Modular Mounting Clamp

**Rail Mount (PC010)** shown on Rotating Shaft (Fits Permobil Unitrack and Pride T-Slotted armrests)

**Tube Mount (PC011)** shown below (Fits Invacare Ø22 mm, 7/8" round-tube armrests)

## 1 SELECT JOYSTICK MOUNTING ADAPTERS



**Quick-Release Lever (FA171)**  
Easily adjust and remove Tri-Lock™ Shaft

### TRI-LOCK™ ROTATING SHAFT (PC024)

Navigate power chairs with smooth precision and prevent accidental rotation from driving position

## 3 SELECT MIDLINE ARM

**Midline Arm**  
Straight arm model shown

### ADDITIONAL KIT OPTIONS

Made of lightweight and corrosion-resistant black anodized aluminum Midline Arms position joysticks where the user needs them to be: above or below arm rests, in driving position or over the shoulder on chin or head-controlled systems. Position joysticks over lap trays, above or below armrests, in driving position or over the shoulder on chin-controlled systems. Bend or cut to size, and support up to 13.6kg (30lbs) of communication device and joystick weight. Available in straight, J- or L-Shapes.




## TODD

"As a man with C4 quadriplegia, my world is in millimeters. My Bodypoint Midline Joystick gives me independence allowing me to drive with my chin – this is critical to my success! I've had the Midline now for more than a decade and it has not failed. Period."




### Tri-Lock™ Rotating Shaft Positions


Shown with a Straight Midline Arm Bracket.




0° Horizontal driving position



90° Allows user to transfer



90° Press to disengage



240° Moved for table access




**3" U-Shaped**  
PC101A Fits stem size: Ø3/16"  
PC101B Fits stem size: Ø1/4"



**4" U-Shaped**  
PC102A Fits stem size: Ø3/16"  
PC102B Fits stem size: Ø1/4"



**Ball knob**  
PC103A Fits stem size: Ø3/16"  
PC103B Fits stem size: Ø1/4"



**Straight**  
PC104A Fits stem size: Ø3/16"  
PC104B Fits stem size: Ø1/4"



**Rubber Dome**  
PC107A Fits stem size: Ø3/16"  
PC107B Fits stem size: Ø1/4"



**2" Rubber Mushroom Head**  
PC108AS Fits stem size: Ø3/16"  
PC108BS Fits stem size: Ø1/4"

### Size Guide

3/16" Stem size:  
Fits most Invacare joysticks  
1/4" Stem size:  
Fits Permobil, Pride, Quantum and Quickie joystick models



**2-1/2" Rubber Mushroom Head**  
PC108AL Fits stem size: Ø3/16"  
PC108BL Fits stem size: Ø1/4"



**Aluminum Top**  
PC109A-BLUE Fits stem size: Ø3/16"  
PC109B-BLUE Fits stem size: Ø1/4"