

# PivotFit™ MD Series Shoulder Harness

Promotes upright posture for improved stability, breathing, and head control. Patented, body-contouring design distributes pressure across a greater surface area than a traditional H-style harness, with balanced retraction from the shoulders down to the lower ribs. Swivel Buckle lets pads pivot to accommodate asymmetry.

Unique shape contours to the shoulders and ribcage for a more supportive, comfortable fit.



MD (multi-directional) Pull Straps can be configured as front-pull or rear-pull at the time of the fitting to meet user needs.



High-strength Swivel Buckle allows pads to pivot into the best placement for the user's body shape, easily accommodating asymmetrical postures and equalizing tension for an ideal fit.



Easy-to-adjust sternum strap is placed low to reduce risk of unintended airway restriction.

## Style Options



### Dynamic Style

Patented, dual-zone stretch enhances comfort by allowing moderate shoulder movement without reducing support at the chest. Ideal for helping to manage spastic tone.



### Standard Style

Fixed, non-stretch style reduces movement and has an easy-to-clean CORDURA® nylon covering. Effective for severe kyphosis, spinal curvature or trunk rotation.

## Part Numbers

Includes four end-fittings. Choose between high-strength Plastic Flat-Mounts and Cinch-Mounts, then add suffix -B5 or -B6 to the end of the part number.

Size	PivotFit™ Dynamic	PivotFit™ Standard
XXS	SH295XXS	SH290XXS
XS	SH295XS	SH290XS
S	SH295S	SH290S
M	SH295M	SH290M
L	SH295L	SH290L
XL	SH295XL	SH290XL

## End-Fittings

### Plastic Flat-Mount™



- B5

### Plastic Cinch-Mount™



- B6

