

Shoulder Harness Strap Guides

Improve positioning and comfort of chest and shoulder harnesses, and use harnesses where you couldn't before — on low or contoured backs. Simple, easy-to-install solution provides maximum adjustability.

Elevate position of top straps to the ideal pull angle for more even pressure distribution over shoulder. Optimize fit for each user.

Super strong spring steel absorbs shock, giving users greater comfort and preventing backrest damage.

Thin design fits into tight spaces.

Integrated holes allow easy harness attachment without extra drilling.

Long mounting slots accommodate a variety of backrests and increase adjustability.

Narrow slot width holds 2.5cm (1") straps in precise position - no sliding off shoulder.

Adjust height, horizontal position, and angle of strap attachment for each side independently.

Compatible with front- and rear-pull strap styles.

MOUNTING:

Attach Strap Guides directly to backrest, as shown at left, or use optional Mounting Adapter plates, shown below.

Strap Guide Options

Available in two lengths to accommodate different backrest sizes and shapes. Sold in pairs with complete M6 fastener kit* and spacers (to accommodate wrapped upholstery or edge trim).



HW100-200-2
Long (200mm/8") Strap Guides



HW100-150-2
Short (150mm/6") Strap Guides

Mounting Adapters (optional)

Strap Guide Mounting Adapters allow you to utilize existing headrest holes, minimizing the need for drilling without limiting options for Strap Guide position. They also make fine-tuning after attachment a snap.

Sold separately, with M6 fasteners and nut plates for harness attachment.



HW105L-1
Large Mounting Adapter



HW105S-1
Small Mounting Adapter

* Standard fasteners compatible with most backrests. For Permobil Corpus chairs, order the HW101P Fastener Pack, sold separately.