



**PD401**  
**SUB-ASIS PAD, 4-5/8" (116mm)**  
**FITS SIZE M, L VELCRO**  
**COMPATIBLE HIP BELTS**

**PD402**  
**SUB-ASIS PAD, 4" (100mm)**  
**FITS SIZE XS, S VELCRO**  
**COMPATIBLE HIP BELTS**

*note*

*DIMENSION ARE IN MILLIMETERS  
 [ALTERNATE UNITS] ARE IN INCHES*

description

**SUB-ASIS PAD SIZES**

|                    |                 |          |          |
|--------------------|-----------------|----------|----------|
| drawing no.        | date originated | drawn by | rev.     |
| <b>D07-0103-01</b> | 1/3/2007        | M. KOSH  | <b>B</b> |