



Instruction: 385 CS Padded Seat Option

Approval: T. Boudreau

Document Number:	SWI-07-031
Revision No.:	4
Revision Date:	10-Jul-14
Written By:	A. Myers

385 CS Padded Seat – Peri Small Hole



Used On: 385 CS Models; **Inception, May 11 2006**
18" & 22" Product Code: SEAT-PD-PERI-SM-CS385
Opening Length: 10.25", Front Width: 4", Rear Width: 6"

385 CS Padded Seat – Peri Large Hole



Used On: 385 CS Models; **Inception, May 11 2006**
Product Code: SEAT-PD-PERI-LG-CS385
Opening Length: 13", Front Width: 4.", Rear Width: 6.5"
Available for 22" width chair only

385 CS Padded Seat – Standard Hole



Used On: 385 CS Models; **Inception, May 11 2006**
18" & 22" Product Code: SEAT-PD-CS385
22" chair: Opening Length: 14", Front Width: 4.5", Rear Width: 6"
18" chair: Opening Length: 13.5", Front Width: 3.5", Rear Width: 5"

385 CS Padded Seat – Transport Solid Seat



Used On: 385 CS Models; **Inception, May 11 2006**
18" & 22" Product Code: SEAT-PD-TS-CS385
Solid Seat 22"