


Encore

Rehab Wheelchair

REHAB
BRODA[®]



Back canes fold down for easy vehicle loading.



Height-adjustable and removable armrests accommodate a variety of users.



Front pivot seat tilt maintains a proper position for self-propulsion.



Open seating system allows for almost any padding configuration.

Patient, Facility, and Caregiver Benefits

- Universal seating system gives flexibility to add/adjust accessories.
- Increases quality of life and engagement with surroundings.
- Allows for easy self-mobility and safer patient maneuverability.
- Reduces risk of injury and falls with supportive stable positioning.

*The Encore of Comfort,
Functionality, and Adaptability.*

SPECIFICATIONS

Tilt - Infinitely Adjustable	up to 33° (combined with rocking)
Weight Capacity	up to 250 lbs. / 158 kg
Seat Widths Available	16", 18", 20", 22" (40cm, 45cm, 50cm, 55cm)
Seat Depths Available	16", 17", 18", 19", 20" (41cm, 43cm, 46cm, 48cm, 51cm)
Back Canes Available	16", 18", 21", 24" (41cm, 46cm, 53cm, 61cm) Angle adjustable to 5°, 10°, 15°, 20°
Floor-to-Seat Heights Available	Adjustable 14" - 19" (36cm - 48cm)
Wheels - Casters	5" rear total lock (12.5cm) 5" rear aluminum non-locking
Wheels - Mag Wheels	20" or 22" front mount (50 cm or 55 cm)
Lifetime Warranty Frame	10 years on 16 Gauge tubular steel frame
Medicare HCPCS Code	E1161, plus additional coded accessories

OTHER PRODUCTS IN THE ENCORE FAMILY

Encore Pedal Wheelchair



Shown with short back, headrest, lateral supports, and front mag wheels.



Shown with long back and shoulder bolsters.

FRAME OPTIONS

