

PO: _____	Date: _____	Customer #: _____
Customer Name: _____	Ship To Name: _____	
Address: _____	Address: _____	
Phone: _____	Phone: _____	
Contact: _____	Contact: _____	
Email: _____	Tag: _____	



Standard Features	
• Comfort Tension Seating®	• 350lb Weight Capacity
• Up to 40° Infinitely Adjustable Seat Tilt	• Trendelenburg Capable
• Up to 90° Infinitely Adjustable Back Recline	• Elevating and Articulating Legrest
• 2 x 6" Rear Multi-Directional Locking Casters	• 2 x 5" Total Lock Front Casters
• Height Adjustable, Swing-Away Removable Arms	
• Height and Width Adjustable Upper Laterals (Wings)	

Synthesis Tilt Recliner - Model: V4 MSRP \$3,674

Seat and Back Options			
Seat Width:	<input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20" N/C	Seat Depth:	<input type="radio"/> 17" Standard
	<input type="radio"/> 22" \$418		<input type="radio"/> 18.5" N/C
	<input type="radio"/> 24" <input type="radio"/> 26" \$946		<input type="radio"/> 20" \$100
Seat Height:	<input type="radio"/> 19.5" \$180	Back Height:	<input type="radio"/> 32" Standard
	<input type="radio"/> 20.5" Standard		<input type="radio"/> 33.5" \$100
	<input type="radio"/> 22" <input type="radio"/> 23.5" \$180		<input type="radio"/> 35" \$100

Front Options	Wheel Options	Rear Options
<input type="radio"/> 5" Total Lock Front Casters Standard	<input type="radio"/> 6" Multi-Directional Locking Rear Casters Standard	
<input type="radio"/> 5" Aluminum Non-Locking Front Casters \$190	<input type="radio"/> 14" Rear Mags <input type="radio"/> 14" Mid Mags \$332	
<i>(available with rear mags only)</i>	<input type="radio"/> 20" Rear Mags \$332	
<input type="radio"/> 5" Front Locking Twin Wheel Casters \$190	<input type="radio"/> Wheel Lock Extension Handles (pair) \$64	

Elevating Footrest Options (Select One)	
<input type="radio"/> Adjustable Elevating Split Footrest (not available in 16" seat width)	Standard
Choose Number of Straps (footrest length range)	<input type="radio"/> 2 (12"-19") <input type="radio"/> 3 (14"-21") <input type="radio"/> 4 (16"-23") <input type="radio"/> 5 (18"-25")
<input type="radio"/> Adjustable Elevating Single Footrest	N/C
Choose Number of Straps (footrest length range)	<input type="radio"/> 2 (12"-19") <input type="radio"/> 3 (14"-21") <input type="radio"/> 4 (16"-23") <input type="radio"/> 5 (18"-25")
<input type="radio"/> ABS Elevating Leg & Footrest (Includes calf and sole pads) (16" - 24" Widths)	\$344

Armrest Options			
Armrest	<input type="radio"/> Soft Armrest Standard	Armtop	<input type="radio"/> Height Adjustable Arms Standard
Options:	<input type="radio"/> Long Armrest \$42	Style:	<input type="radio"/> Height and Depth Adjustable (pair) \$138
	<input type="radio"/> Padded Armrests (pair)(black) \$80	Armtop	<input checked="" type="radio"/> Seat Mounted Arms Standard
	<input type="radio"/> Arm Trough (each) \$136	Mount:	

Base Frame Finish Color (Select one)			
<input type="radio"/> Carmine Red	<input type="radio"/> Ocean Blue	<input type="radio"/> Midnight Black	N/C
<input type="radio"/> Satin Nickel	<input type="radio"/> Tungsten Sand	<input type="radio"/> Bronze	N/C

Upholstery Color Options (Select one)

- * Standard padding configuration: Headrest, exterior of shoulder bolsters, upper and lower side pads & sole pads Come in Accent Color. Seat, back and calf pads, and interior of shoulder bolsters come standard in Savoy Black Vinyl.
- ** Indigo, Espresso & Graphite are available on seat, back & calf pad & most accessories. See below under padding for selection.
- *** In honor of our Veterans Service, Camo options are now available.



Standard Padding*: Shown with Geyser Upholstery & Tungsten Sand Frame Color

Accent Padding:** Shown with Indigo Upholstery & Black Frame Color *Chair Shown with ABS footrest

Upholstery Accent Color* (Select One)

- | | |
|--|--|
| Vinyl: | Polyurethane: |
| <input type="radio"/> Indigo** <input type="radio"/> Interlace Espresso** <input type="radio"/> Graphite** | <input type="radio"/> Portobello <input type="radio"/> Geyser <input type="radio"/> Fossil <input type="radio"/> Forge |
| <input type="radio"/> Black <input type="radio"/> Smoke House <input type="radio"/> Chocolate | Camo Nylon:*** (with performance finish) |
| | <input type="radio"/> Digital Camo Green <input type="radio"/> Retro Camo Jungle \$130 |

Padding - Choose Option 1 OR 2

- Option 1 - BRODA Padding** N/C

Standard BRODA Padding:

Includes: **In Accent Color:** Headrest, Exterior Shoulder Bolsters, Upper & Lower Side Pads & Sole Pads*.

In Black Vinyl: Seat & Back Pads with Covers, Calf Pad & Interior Should Bolsters

Seat, Back & Calf	<input type="radio"/> Black Vinyl	Standard
Pad Cover Options:	<input type="radio"/> Accent Color*	N/C
	*** <i>(For Indigo, Espresso or Graphite only)</i>	
	<input type="radio"/> Full vinyl Wrap Seat & Back Pads	\$100
	<i>(Select Black Vinyl OR Accent Color Vinyl Above)</i>	
	<input type="radio"/> Calf Pad Color	N/C
	<i>(Select Black Vinyl OR Accent Color Vinyl Above)</i>	
	<input type="radio"/> Grey Terry <input type="radio"/> White Terry	N/C
HSP Padding Package:	<input type="radio"/> Includes Shoulder Bolsters Long, Thigh Belt & Footbox Long	\$544
	T-Belt Color: <input type="radio"/> Black Vinyl <input type="radio"/> Accent Color Vinyl	N/C
	Footbox Long Color: <input type="radio"/> Black Vinyl <input type="radio"/> Accent Color Vinyl	N/C
	<input type="radio"/> Optional Elevating Fixed Footrest	N/C

- Option 2 - BRODA APP (Additional Positioning Package)**

BRODA APP Padding Package Includes: Seat Cushion, Back and Head Support, Foot & Calf Support, Upper & Lower Side Pads

Seat Width:	<input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20"	Seat Depth:	<input type="radio"/> 17" <input type="radio"/> 19" <input type="radio"/> 21"	\$522
	<input type="radio"/> 22" <input type="radio"/> 24"		<input type="radio"/> 17" <input type="radio"/> 19" <input type="radio"/> 21"	\$522
	<input type="radio"/> 26"		<input type="radio"/> 17" <input type="radio"/> 19" <input type="radio"/> 21"	\$522

- Acute Care Padding Package:**

Includes: Full Vinyl Covers, ABS Footrest with Calf and Sole Pads *(Factory installed)(skip footrest and padding package section above when selecting this option)(16"-24" widths)* \$476

Common Accessories					
Seat Belts:	<input type="radio"/> S-BELT	Seat Belt	<input type="radio"/> Center Buckle	<input type="radio"/> Side Buckle	\$52
	<input type="radio"/> S-BELT-A	Seat Belt Auto Buckle	<input type="radio"/> Center Buckle	<input type="radio"/> Side Buckle	\$52
	<input type="radio"/> S-BELT-PD	Seat Belt Padded			\$76
	<input type="radio"/> S-BELT-PD-A	Seat Belt Padded Auto Buckle			\$76
Thigh Belts: <i>(not compatible with BRODA APP padding option)</i>	<input type="radio"/> T-BELT	Padded Thigh Belt			\$122
	<input type="radio"/> T-BELT-A	Padded Thigh Belt Auto Buckle			\$126
	<input type="radio"/> T-BELT-NEO	Padded Thigh Belt Neoprene			\$122
	<input type="radio"/> T-BELT-NEO-A	Padded Thigh Belt Neoprene Auto Buckle			\$126

Common Accessories (Continued)					
Trays:	<input type="radio"/> TRAY-ABS	Tray ABS with Strap	<input type="radio"/> Sewn Tray Strap	<input type="radio"/> Grey <input type="radio"/> Black	\$312
	<input type="radio"/> TRAY-ABS-OVR	Tray ABS with Strap (24" - 26" width)		<input type="radio"/> Black Only	\$448
	<input type="radio"/> TRAY-CLEAR	Tray Clear with Strap	<input type="radio"/> Sewn Tray Strap		\$322
	<input type="radio"/> TRAY-HALF	Half Tray (side specific)(includes upper side pad)			\$428
	<input type="radio"/> TRAY-COVER	Tray Cover (padded)(for TRAY-ABS)			\$158
Stabilizers:	<input type="radio"/> STABIL-FR	<input type="radio"/> Front or <input type="radio"/> Rear Mounted Stabilizers (select one)(pair)		\$90	
	<input type="radio"/> STABIL-LAT	Lateral Frame Stabilizers		\$158	
Other:	<input type="radio"/> IV-POLE	IV Pole & Holder (includes mounting bracket) - <input type="radio"/> Front or <input type="radio"/> Rear Mount		\$100	
	<input type="radio"/> O2-HOLDER	O2 Holder (includes mounting bracket) - <input type="radio"/> Front or <input type="radio"/> Rear Mount		\$168	

Common Padding Upgrades				
Headrests:	<input type="radio"/> NR-C	Headrest Contoured Pad		\$106
	<input type="radio"/> NR-CV	Headrest Visco Foam Contoured Pad	<input type="radio"/> Omit Standard	\$106
	<input type="radio"/> NS	Neck Support	Headrest Less \$64	\$106
Positioning:	<input type="radio"/> SR	Side Reduction Pads (pair)		\$100
	<input type="radio"/> CC-LTS	Comfort Company Lateral Supports		\$472
	<input type="radio"/> LTS	Lateral Torso Supports (FACTORY INSTALLED RECOMMENDED)		\$580
	<input type="radio"/> SB-LONG-V4	Shoulder Bolster Long (pair) - Omits standard short		\$286
Footrest Pads:	<input type="radio"/> FB-SH	Footbox - Extended Pad	<input type="radio"/> Omit Calf and Sole	\$300
	<input type="radio"/> FB-LONG	Footbox- Full Extended Pad	Pads Less \$178	\$344

Additional Options: _____

Delivery Instructions or Special Requests: _____

Total: \$
