

# The Innovative Reclining Manual Wheelchair For Kids, Teens & Young Adults





SPECIFICATIONS	ZREC1200	ZREC1400	ZREC1600
FRAME COLORS	ORANGE	RED	BLUE
SEAT WIDTH	12"	14"	16"
SEAT DEPTH	12"	14"	16"
OVERALL OPEN WIDTH	20"	20"	24"
OVERALL FOLDED WIDTH	12"	12"	12"
BACK HEIGHT INCL HEADREST	25"	25"	31"
BACK HEIGHT W.O HEADREST	17.5"	17.5"	19.5"
OVERALL HEIGHT	40.75"	40.75"	49.5"
SEAT TO FLOOR HEIGHT	15" / 16"	15" / 16"	19.5" / 20.5"
OVERALL LEG W/ RIGGINGS	38.5"	38.5"	52.75"
BACKREST RECLINES	90° - 180°	90° - 180°	90° - 180°
NET WEIGHT	44 LBS	44 LBS	56.75 LBS
CARTON SHIPPING WEIGHT	50 LBS	50 LBS	62 LBS
WEIGHT CAPACITY	150 LBS	200 LBS	300 LBS



Lightweight foam cushion provides comfort and sitting tolerance



**ACCESSORY**

Reclining back with adjustments up to 180 degrees

**ZGCS12** for use with **ZGREC1200**  
**ZGCS14P** for use with **ZGREC1400**  
**ZGCS16** for use with **ZGREC1600**



Expertly designed to provide a whole new found level of freedom, independence and comfort by combining:

- ✓ A Range Of Positioning Options
- ✓ Reclining Back With Infinite Adjustments Up To 180°
- ✓ Dual Axle For Seat-To-Floor Height Adjustability
- ✓ Anti-tippers provide stability and prevent tipping.
- ✓ Colored Steel Frame, Lightweight & Foldable Design

It's time for your child to enjoy the gift of easier, safer and more comfortable movement, indoors and outdoors.

# Dream Big Roam Free

This Time, with

**ziggo**Pro



Powered by Circle Specialty's 25+ years of experience in revolutionary mobility and rehab equipment tailored exclusively to the needs of kids and teens, you have every good reason to make **ziggo**Pro a part of your family.