

# Seatar


## BathMobile

Designed for you  
commode & shower chair



seat turn 180°

adjustable footrest



### features



Weighs 10 kg.

The seat opening  
can be rotated  
180 degrees.

Drive wheels can  
be assembled.

Long armrests  
that assist in  
seating and rising.

Adjustable  
footrest.

Detachable  
bedpan rails.

The chair is made  
with advanced plastic  
injection technology.

All metal parts are  
made of stainless steel  
(including bearings).





Seatarra

# BathMobile

Designed for you

**BathMobile** is a breakthrough product in a conservative world.

**BathMobile** is a modular product that can be detached and assembled without tools, making it easy to move around.

**BathMobile** allows you the comfort you are used to elsewhere.

**BathMobile** allows you the freedom to change and adjust; swivel seat, preparation for drive wheels, adjustable product height, detachable bedpan rails.



**Clarke**  
Health Care Products

Distributed by:  
**Clarke Health Care Products**  
7830 Steubenville Pike • Oakdale, PA 15071