

Chrome Two-Step Step Stool



MODEL: T-6842

- Durable chrome finish
- Extra-wide and deep tops
- Reinforced base and top
- First step 10" high x 12.25" x 17.25"
- 350 lbs. load capacity under normal use

**Shipping weight: 30 lbs.*

Box dimensions: 29" x 19" x 19"

** This item could be dimensionally weighted when shipped by UPS® or FedEx® Ground*

SPECIFICATIONS

PLATFORM: 14.25" L X 11.25" W X 9"H



California residents please be advised of the following, pursuant to Proposition 65:



WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity.

1

Chrome Two-Step Step Stool with Handrail



MODEL: T-6850

- Angled handrail 43" at highest point
- Durable chrome finish
- Extra-wide and deep tops
- Reinforced base and top
- First step 10" high x 12.25" x 17.25"
- 350 lbs. load capacity under normal use

NOTE: Handrail can be attached to either side of the stool for left hand or right hand use.

**Shipping weight: 38 lbs., in 2 boxes*

Box dimensions: 29" x 19" x 19" and 43" x 27" x 2"

SPECIFICATIONS

27" W X 17.25" D X 17.25"H



California residents please be advised of the following, pursuant to Proposition 65:



WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity.

1