

## Large Top Bariatric Step Stool



### MODEL: T-6242

- Durable chrome finish
- Extra wide and deep top
- Welded reinforced "X" brace
- 600 lbs. load capacity under normal use


*Shipping weight: 14 lbs.*

### SPECIFICATIONS

PLATFORM: 17.25" L X 12.25" W X 9"H



California residents please be advised of the following, pursuant to Proposition 65:

 **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. **1**

## Large Top Bariatric Step Stool with Handrail



### MODEL: T-6250

- 1" diameter handrail 36" high
- Durable chrome finish
- Extra wide and deep top
- Welded reinforced "X" brace
- 600 lbs. load capacity under normal use


*Shipping weight: 18 lbs.*

### SPECIFICATIONS

PLATFORM: 17.25" L X 12.25" W X 9"H



**California residents please be advised of the following, pursuant to Proposition 65:**

 **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. **1**