

# Accessories for Toilet Support System

## Armrest and Footrest Options

**MODEL #4001** - U-Shape Mounting Base only for Toilet Support System  
18½" W x 3" D x 10" H

**MODEL #4921** - Armrest only, no U-Shape Mounting Base

**MODEL #4922** - Footrest only, no U-Shape Mounting Base

**MODEL #4925** - Armrest and Footrest with brackets, no U-Shape Mounting Base

**MODEL #4930** - Armrest and U-Shape Mounting Base only

Armrest Width Outside: 19"  
Armrest Width Inside: 15½"  
Height Adjustment: 20" - 31"  
Length Adjustment: 19" - 22"  
Footrest Dimensions: 17" W x 7" H x 7" D  
Leg Drop: 17" to floor



## Plastic Reducer Ring

Fits on any toilet seat instantly. Sanitary, seamless, easy to clean. Build-in splashguard.

**MODEL #4950**

**Hole:** 6" W x 7" L

**Outside:**  
10" W x 13" L

**Abductor:**  
3½" W x 2" H



## Padded Reducer Ring

For those who require extended toileting time. Includes replaceable abductor splashguard. Sanitary, seamless, easy to clean. Fits on any toilet seat instantly, 1½" thick.

**MODEL #4960**

**Hole:** 7" W x 8" L

**Outside:**  
14" W x 15½" L

**Abductor:**  
3½" W x 4½" H

**MODEL #4965**

**Hole:** 5" W x 6" L

**Outside:**  
14" W x 15½" L

**Abductor:**  
3½" W x 3½" H



## Soft-Flex™ Splashguard

For adults and larger children. Slides off easily after use. Fits open or closed front toilet seat. Cleans easily.

**MODEL #4975**

5" W x 5½" L

**MODEL #4976**

Replacement Straps  
2 ea. for Splashguard