

COMFORT FOOT SINGLE



IST - Inferior Support Thickness
L - Posterior Support Length

PRODUCT FEATURES

- Unique heel cut-out eliminates contact with heel, prevents pressure points, and reduces risk of skin breakdown.
- Optional gel pack in support base reduces pressure, shear, and heat build-up.
- Medium density, high resiliency foam construction provides maximum comfort.
- Non-skid fabric helps prevent unwanted foot migration.
- Tapered bottom accommodates flip up foot supports.
- Fits on both left and right foot.

HARDWARE

Clamp-On Connection: Easy on, easy off connection gives users a versatile attachment.

Bolt-On Connection: Bolt-on connection works with all foot plates. Recommended for users with tone.



COMFORT FOOT SINGLE

Style	W x D x L	Part Number	Gel Part Number
Standard	7 x 10 x 8	FTS	FTSG
Tall	7 x 10 x 12	FTST	FTSTG
Bariatric	10 x 10 x 8	FTSB	FTSBG

INFERIOR SUPPORT PAD (BASE ONLY)

Style	W x D	Part Number	Gel Part Number
Standard	7 x 10	FTS-BTM	FTSG-BTM
Bariatric	10 x 10	FTSB-BTM	FTSBG-BTM

Style	Inferior Support Thickness (IST)
Standard	2"
W/Gel	2.5"
Style	Posterior Support Thickness
Standard	1.75"
W/Gel	1.75"

COMFORT FOOT SINGLE ELEVATING

Style	W x D x L	Part Number
Standard	7 x 10 x 8	FTSE
Tall	7 x 10 x 12	FTSTE

POSTERIOR SUPPORT PAD (TOP PORTION ONLY)

Style	W x D x L	Part Number
Standard Top Portion Only	7 x 10 x 8	FTS-TOP
Tall Top Portion Only	7 x 10 x 12	FTST-TOP
Bariatric Top Portion Only	10 x 10 x 8	FTSB-TOP
Elevating Top Portion Only	7 x 10 x 8	FTSE-TOP
Tall Elevating Top Portion Only	7 x 10 x 12	FTSTE-TOP

Reimbursement for items billed using K0108 is not guaranteed, as coverage is at the carrier's discretion.

Comes standard with Comfort-Tek. See Pg. 2 for details.