

CURVE®

E2601 / E2602
GENERAL USE

PRODUCT FEATURES



1. Mild anti-thrust shelf promotes postural stability and helps prevent forward migration.
2. Medial thigh separator promotes postural alignment.
3. High density, molded foam increases surface area reducing potential peak pressure areas.
4. Non-skid flat bottom fabric offers hook and loop attachment.

DIMENSIONS

Seat Thickness	3"
Inferior Thigh Support Thickness	2.25"
Medial Thigh Support Thickness	1"
Lateral Thigh Support Height	1"
Pelvic Contour Height	1"

WEIGHT CAPACITY

12" Wide & Under	100 lbs
14" Wide	150 lbs
16" Wide	250 lbs
18" Wide	300 lbs
20" Wide	350 lbs
22" Wide & Wider	400 lbs

PEDIATRIC SIZES

	10" W	12" W	14" W
10" D	10 x 10	12 x 10	14 x 10
12" D	10 x 12	12 x 12	14 x 12
14" D	10 x 14	12 x 14	14 x 14
16" D		12 x 16	14 x 16
18" D			14 x 18

CLINICIAN TIP
Top general use choice for pressure redistribution and stability

STANDARD SIZES

Additional Sizes Available

	16" W	18" W	20" W	22" W	24" W
14" D	16 x 14				
16" D	16 x 16	18 x 16	20 x 16	22 x 16	24 x 16
18" D	16 x 18	18 x 18	20 x 18	22 x 18	24 x 18
20" D	16 x 20	18 x 20	20 x 20	22 x 20	24 x 20

CURVE PART NUMBERS


Cover Style	Stretch-Air™	Comfort-Tek™
Standard	CU-SV-WWDD	CU-FV-WWDD

WWDD to be replaced / WW = width & DD = depth in inches

ELEMENTS® WITH GEL

E2601 / E2602
GENERAL USE


PRODUCT FEATURES



1. Four compartment inset QuadraGel® is designed to prevent gel migration and redistributes pressure to provide maximum comfort.
2. Non-skid flat bottom fabric offers hook and loop attachments.
3. Single layer foam construction contours to the user's body reducing high pressure areas.

QUADRAGEL®

The QuadraGel® prevents gel migration and redistributes pressure for comfort



WEIGHT CAPACITY

16" Wide	250 lbs
18" Wide	300 lbs
20" Wide	350 lbs
22" Wide & Wider	400 lbs

STANDARD SIZES

	16" W	18" W	20" W	22" W	24" W
16" D	16 x 16	18 x 16	20 x 16	22 x 16	24 x 16
18" D	16 x 18	18 x 18	20 x 18	22 x 18	24 x 18
20" D	16 x 20	18 x 20	20 x 20	22 x 20	24 x 20

Additional Sizes Available

ELEMENTS WITH GEL PART NUMBERS

Seat Thickness	Stretch-Air™	Comfort-Tek™
2" Seat Thickness	46GS-WWDD-B	46G-WWDD-B
3" Seat Thickness	463GS-WWDD-B	463G-WWDD-B

WWDD to be replaced / WW = width & DD = depth in inches