

## REPLACEMENT COVER

**E2619**  
REPLACEMENT COVER

**GLIDE  
WEAR.**

Available with GlideWear®  
Technology. See pg. 25 for details.



### STRETCH-AIR™

Stretch-Air™ is designed to be airy and breathable with superior stretch. It provides heat dissipation for user comfort while the multi-directional stretch contributes to pressure relief by allowing the user to gain the full benefits of the cushion design.

### COMFORT-TEK™

Comfort-Tek™ is smooth and strong yet extremely soft for comfort. The multidirectional stretch provides pressure relief by allowing the fabric to better conform to the individual's unique shape. Designed for infection control and as an incontinence barrier, Comfort-Tek™ is easily cleaned anytime with common disinfectants and is latex-free.

### STRETCH-AIR™



### COMFORT-TEK™



REPLACEMENT COVERS ARE AVAILABLE FOR  
SEATING SURFACES AND BACK SUPPORTS

## INCONTINENT PROOF LINER

**E2619**  
REPLACEMENT COVER



### PRODUCT FEATURES

- Provides extra incontinence protection for seating surfaces.
- The liner may be used as an extra cover when the outer cover is being washed.
- Cushion style and size must be specified at time of order.

### INCONTINENT PROOF PART NUMBERS

BINC-WWDD-CV

WWDD to be replaced / WW = width & DD = depth in inches