


# ELEMENTS® with GEL

HCPCS CODING: E2601/E2602  
GENERAL USE



## PRODUCT FEATURES




**NEW:** Elements with Gel now comes standard with non-skid bottom and hook/loop connection.



### QUADRAGEL PACK

The QuadraGel pack regulates body temperature and redistributes pressure for comfort.



## COVER OPTIONS



**STRETCH-AIR FABRIC:** designed for air transmission and heat dissipation. This breathable laminate combines a high-stretch polyester top layer with a waterproof, breathable vapor barrier.



**COMFORT-TEK FABRIC:** the fluid proof exterior has supple leathery properties that provide an easily cleaned surface and protection for the inner cushion. Comfort-Tek has multi-dimensional stretch and helps alleviate pressure on bony prominences.

### ELEMENTS WITH GEL

2" THICK

stretch-air part #	comfort-tek part #	w x l	thickness	wt. limit
<b>46GS-1616-B</b>	<b>46G-1616-B</b>	16 x 16	2"	250 lbs
<b>46GS-1618-B</b>	<b>46G-1618-B</b>	16 x 18	2"	250 lbs
<b>46GS-1620-B</b>	<b>46G-1620-B</b>	16 x 20	2"	250 lbs
<b>46GS-1816-B</b>	<b>46G-1816-B</b>	18 x 16	2"	300 lbs
<b>46GS-1818-B</b>	<b>46G-1818-B</b>	18 x 18	2"	300 lbs
<b>46GS-1820-B</b>	<b>46G-1820-B</b>	18 x 20	2"	300 lbs
<b>46GS-2016-B</b>	<b>46G-2016-B</b>	20 x 16	2"	350 lbs
<b>46GS-2018-B</b>	<b>46G-2018-B</b>	20 x 18	2"	350 lbs
<b>46GS-2020-B</b>	<b>46G-2020-B</b>	20 x 20	2"	400 lbs
<b>46GS-2216-B</b>	<b>46G-2216-B</b>	22 x 16	2"	400 lbs
<b>46GS-2218-B</b>	<b>46G-2218-B</b>	22 x 18	2"	400 lbs
<b>46GS-2220-B</b>	<b>46G-2220-B</b>	22 x 20	2"	400 lbs
<b>46GS-2416-B</b>	<b>46G-2416-B</b>	24 x 16	2"	400 lbs
<b>46GS-2418-B</b>	<b>46G-2418-B</b>	24 x 18	2"	400 lbs
<b>46GS-2420-B</b>	<b>46G-2420-B</b>	24 x 20	2"	400 lbs

NOTES

+ CUSHIONS AVAILABLE WITH RIGID INSERTS & EXTRA INCONTINENT PROOF LINER. SEE PAGE 19 FOR DETAILS

### ELEMENTS WITH GEL

3" THICK

stretch-air part #	comfort-tek part #	w x l	thickness	wt. limit
<b>463GS-1616-B</b>	<b>463G-1616-B</b>	16 x 16	3"	250 lbs
<b>463GS-1618-B</b>	<b>463G-1618-B</b>	16 x 18	3"	250 lbs
<b>463GS-1620-B</b>	<b>463G-1620-B</b>	16 x 20	3"	250 lbs
<b>463GS-1816-B</b>	<b>463G-1816-B</b>	18 x 16	3"	300 lbs
<b>463GS-1818-B</b>	<b>463G-1818-B</b>	18 x 18	3"	300 lbs
<b>463GS-1820-B</b>	<b>463G-1820-B</b>	18 x 20	3"	300 lbs
<b>463GS-2016-B</b>	<b>463G-2016-B</b>	20 x 16	3"	350 lbs
<b>463GS-2018-B</b>	<b>463G-2018-B</b>	20 x 18	3"	350 lbs
<b>463GS-2020-B</b>	<b>463G-2020-B</b>	20 x 20	3"	400 lbs
<b>463GS-2216-B</b>	<b>463G-2216-B</b>	22 x 16	3"	400 lbs
<b>463GS-2218-B</b>	<b>463G-2218-B</b>	22 x 18	3"	400 lbs
<b>463GS-2220-B</b>	<b>463G-2220-B</b>	22 x 20	3"	400 lbs
<b>463GS-2416-B</b>	<b>463G-2416-B</b>	24 x 16	3"	400 lbs
<b>463GS-2418-B</b>	<b>463G-2418-B</b>	24 x 18	3"	400 lbs
<b>463GS-2420-B</b>	<b>463G-2420-B</b>	24 x 20	3"	400 lbs

NOTES

+ CUSHIONS AVAILABLE WITH RIGID INSERTS & EXTRA INCONTINENT PROOF LINER. SEE PAGE 19 FOR DETAILS