

# COMFORT ARM<sup>®</sup>



STANDARD ARM



LOW ARM



TALL RAIL ARM



MOLDED ARM


## E2209 ARM SUPPORT

## E1028 ARM SUPPORT HARDWARE

### PRODUCT FEATURES

- Comfort Arm design provides finger support for all users.
- Comfortable straps for securing the arm can be adjusted or removed.
- Provides support for contractures and tone while also helping to prevent finger drop.
- Moldable option can be used to elevate fingers.

### COMFORT ARM DIMENSIONS



D - Total Support Depth  
W1 - Total Support Width  
W2 - Effective Support Width  
T1 - Elbow Well Height  
T2 - Inferior Support Thickness

### EAD HARDWARE

U.S. Patent No. 7,866,613

Allows arm support to be elevated and articulated with a quick release lever.

Arm support can be adjusted for a depth of up to 2.5" and a width of up to 4".

Hardware weighs only 1.2 lbs.



ELEVATING



ARTICULATING



DEPTH-ADJUSTABLE

### INTERLOCK BRACKET

The interlock bracket locks down to 3/4" up to 1" tubing with a two-bolt connection.

Standard with fixed or swing-away hardware option.

### SWING-AWAY HARDWARE

The rotating base system allows the Comfort Arm to be angularly rotated side-to-side and adjusted front to back 7".

### SLIDE ON HARDWARE

Slides arm support over stock wheelchair arm pad.

### TB3 HARDWARE

Arm support with track mount hardware for TB3 style chairs.

### COMFORT ARM SUPPORTS

STANDARD ARM SUPPORTS					
D x W1	W2	T1	T2	Side	Model Number
22" x 6"	4.0"	2.75"	1.5"	Right	2002R-PAD
22" x 6"	4.0"	2.75"	1.5"	Left	2002L-PAD
14" x 6"	4.0"	2.75"	1.5"	Right	2002R-14-PAD
14" x 6"	4.0"	2.75"	1.5"	Left	2002L-14-PAD
17" x 6"	4.0"	2.75"	1.5"	Right	2002R-17-PAD
17" x 6"	4.0"	2.75"	1.5"	Left	2002L-17-PAD

STANDARD ARM MOLDABLE SUPPORTS					
D x W1	W2	T1	T2	Side	Model Number
22" x 6"	4.0"	2.75"	1.5"	Right	2002R-M-PAD
22" x 6"	4.0"	2.75"	1.5"	Left	2002L-M-PAD

TALL RAIL ARM SUPPORTS					
D x W1	W2	T1	T2	Side	Model Number
22" x 7"	4.0"	2.75"	1.5"	Right	2007R-PAD
22" x 7"	4.0"	2.75"	1.5"	Left	2007L-PAD

LOW ARM SUPPORTS					
D x W1	W2	T1	T2	Side	Model Number
22" x 6"	3.0"	n/a	1.25"	Right	2004R-PAD
22" x 6"	3.0"	n/a	1.25"	Left	2004L-PAD
14" x 6"	3.0"	n/a	1.25"	Right	2004R-14-PAD
14" x 6"	3.0"	n/a	1.25"	Left	2004L-14-PAD
17" x 6"	3.0"	n/a	1.25"	Right	2004R-17-PAD
17" x 6"	3.0"	n/a	1.25"	Left	2004L-17-PAD

LOW ARM MOLDABLE SUPPORTS					
D x W1	W2	T1	T2	Side	Model Number
22" x 6"	3.0"	n/a	1.25"	Right	2004R-M-PAD
22" x 6"	3.0"	n/a	1.25"	Left	2004L-M-PAD

TALL RAIL ARM MOLDABLE SUPPORTS					
D x W1	W2	T1	T2	Side	Model Number
22" x 7"	4.0"	2.75"	1.5"	Right	2007R-M-PAD
22" x 7"	4.0"	2.75"	1.5"	Left	2007L-M-PAD

STANDARD ARM FOR TB3					
D x W1	W2	T1	T2	Side	Model Number
22" x 6"	4.0"	2.75"	1.5"	Right	1002R-TM
22" x 6"	4.0"	2.75"	1.5"	Left	1002L-TM
14" x 6"	4.0"	2.75"	1.5"	Right	1002R-14-TM
14" x 6"	4.0"	2.75"	1.5"	Left	1002L-14-TM
17" x 6"	4.0"	2.75"	1.5"	Right	1002R-17-TM
17" x 6"	4.0"	2.75"	1.5"	Left	1002L-17-TM

Moldable Option					
D x W1	W2	T1	T2	Side	Model Number
22" x 6"	4.0"	2.75"	1.5"	Right	1002R-M-TM
22" x 6"	4.0"	2.75"	1.5"	Left	1002L-M-TM

### COMFORT ARM HARDWARE

COMFORT ARM HARDWARE		
Hardware Type	Side	Model Number
Fixed Interlock	Right	1000R-HARD
Fixed Interlock	Left	1000L-HARD
Swing Away	Right	2000R-HARD
Swing Away	Left	2000L-HARD
EAD	Right	3000R-HARD
EAD	Left	3000L-HARD
Slide On	Right	2100R-HARD
Slide On	Left	2100L-HARD
TB3 Track Mount	N/A	ARMTM-P3-HARD

Comes standard with Comfort-Tek. See Pg. 2 for details.

### MOLDED ARM

MOLDED ARM SUPPORTS		
Depth	Side	Model Number
22"	Right	2300R-PAD
22"	Left	2300L-PAD
14"	Right	2300R-14-PAD
14"	Left	2300L-14-PAD

MOLDED ARM FOR EAD HARDWARE		
Depth	Side	Model Number
22"	Right	2300PR-PAD
22"	Left	2300PL-PAD
14"	Right	2300PR-14-PAD
14"	Left	2300PL-14-PAD