

Comfy™ Grip Hand

Gradual finger extension, contracture management, resting comfort

SKIN SAFE

Anti-microbial material for skin breakdown prevention

ADJUSTABLE

Patented reinforced spine bends to wrist and forearm needs

PRACTICAL

Holds hand in a functional C position and helps control wrist drop

CONTRACTURE MANAGEMENT

Frame designed to slide between palm and finger of even the most contracted hands

EXTENSION

Two graduated rolls gradually increase extension of the MCP and PIP joints

SIDE WINGS

Adjustable wings prevent ulnar and radial deviation

COMFORT

Super Soft and removable Headliner cover is machine washable



L3807

L3809

	Wrist to Palmer Creases	Wrist Down Forearm	Total Splint Length
Large	3 ¼"	7"	10 ¼"
Medium	2 ¾"	6 ¼"	9"
Small	2"	5"	7"