

Comfy™ Hand Finger Wrist Orthosis

A therapist favorite for good reason...



Broadcloth Comfy cover

Ambidextrous

RELAXING

An ideal resting splint

SIDE WINGS

Prevents ulnar and radial deviation

POSSIBILITIES

Comes in many variations:

- Comfyprene- Slim line, no cover
- Long Hand
- Long Arm
- Pediatric
- Broadcloth or Terrycloth cover
- 4 straps



Comfyprene options shown

ADJUSTABLE

Patented spine bends to wrist and forearm needs

COMFORT

Super Soft and removable Headliner or Terrycloth covers are machine washable

SKIN SAFE

Anti-microbial and moisture wicking material for skin breakdown prevention



	Wrist to Tip of Middle Finger	Total Length	Comfy
Adult Large	8" and Above	13 ¾"	H-101-AL
Adult	7" – 8"	12 ½"	H-101
Adult Small	6" – 7"	11"	H-101-AS
Pedi. Lg	4 ½" – 6"	9 ½"	PH-101-L
Ped. Md	3" – 4 ½"	8"	PH-101-M
Ped. Sm.	2 ½" – 3 ½"	7"	PH-101-S