

# Comfy™ Hand Thumb Orthosis

Our Therapist-loved Hand Wrist orthosis with Extra Thumb Support

## SKIN SAFE

Anti-microbial and moisture wicking material for skin breakdown prevention



## COMFORT

Super Soft and removable Headliner or Terrycloth covers are machine washable

## THUMB SUPPORT

Unique folding thumb tabs allow for extra thumb control and support.



Shown in Navy Comfypene version

## ADJUSTABLE

Patented spine bends to wrist and forearm needs

## RELAXING

An ideal resting splint

## POSSIBILITIES

Comes in many variations:

- Comfyprene- Slim line, no cover
- Long Hand
- Long Arm
- Pediatric
- Broadcloth or Terrycloth cover
- 4 straps

## SIDE WINGS

Prevents ulnar and radial deviation

Size	Wrist to Tip of Middle Finger	Total Length	COMFY #
Adult Large	8" and Above	13 ¾"	HT-101-AL
Adult	7" – 8"	12 ½"	HT-101
Adult Small	6" – 7"	11"	HT-101-AS
Pediatric Large	4 ½" – 6"	9 ½"	PHT-101-L
Pediatric Medium	3" – 4 ½"	8"	PHT-101-M
Pediatric Small	2 ½" – 3 ½"	7"	PHT-101-S

L3807  
or  
L3809

Ambidextrous