

Comfy™ Adjustable Cone Hand Orthosis

Patented adjustable splint to treat strong flexion synergy

- Lateral forearm support
 - Therapist Loved
- Allows for wrist movement



Shown with included removable I finger separator

SKIN SAFE

Anti-microbial material for skin breakdown prevention

COMFORT

Super Soft and removable Headliner cover is machine washable

MOVABLE

Adjustable hinge at wrist helps increase wrist extension or flexion as needed



FINGER CONE

Worn laterally with fingers positioned around cone. Adjusts to contractures



SIDE WINGS

Adjustable wings prevent ulnar and radial deviation

ADJUSTABLE

Patented spine bends along with movable wrist joint

Ambidextrous

High L3915 (or L3916) Reimbursement

	Wrist to Tip to Middle Finger	Total Length	COMFY
Adult	7" – 8"	12 ½"	ACH-101
Pediatric Lg	4 ½" – 6"	9 ½"	PACH-101-L
Pedi. Med	3" – 4 ½"	8"	PACH-101-M
Pediatric Sm.	2 ½" – 3 ½"	7"	PACH-101-S