

Comfy™ Multi Podus Contracture Resting and Ambulating Boots

SKIN SAFE

Treated with Antimicrobial surface protector to shield against mold and mildew, stains, bacteria fungus and algae

HEEL FLOAT

Boots position heel in floating position to prevent pressure areas

AMBULATE

Ambulation version (A-Boot) comes with removable, soft, non-slip sole

Sample Indications:

- ✓ plantar flexion
- ✓ foot drop
- ✓ wound care
- ✓ ankle edema
- ✓ impaired mobility
- ✓ prevention/treatment of heel/ankle/foot pressure sores



L4396 or
L4397

COMFORT

Super Soft fleece lined Headliner cover is removable and machine washable

ANTI-ROTATION

Rotator bar controls hip and leg rotation

EXPANDABLE

Comes with extender kit to accommodate for wider or edematous feet



ADJUSTABLE TOE POST

Extend toe post to accommodate all foot sizes- universal sizing

C-Boot

Comfy Multi-Podus Boot
(without ambulation)

A-Boot

Comfy Multi-Boot
with Ambulation pad