

Comfy™ Boot Strap

Comfy™ Ambulation Boot Strap

SKIN SAFE

Treated with Antimicrobial surface protector to shield against mold and mildew, stains, bacteria fungus and algae



COMFORT

Super Soft fleece lined Headliner cover is removable and machine washable

HEEL FLOAT

Boots position heel in floating position to prevent pressure areas



STRAPS

Dorsiflexion straps to prevent plantar flexion

ADJUSTABLE TOE POST

Extend toe post to accommodate all foot sizes



EXPANDABLE

Comes with extender kit to accommodate for wider or edematous feet

AMBULATE

Ambulation version comes with removable, soft, non-slip sole

SUPPORT

Ankle strap secures to prevent ankle flexion

ANTI-ROTATION

Rotator bar controls hip and leg rotation

Reduce readmission risk

L4396

	Adult	Pediatric
Foot Length	8" to 13"	5" to 8"
Calf Circumference	Up to 16"*	Up to 11" *
Heel to Top of Back	11 ½"	7"