

Comfy™ Adjustable Cone Hand Orthosis

Patented adjustable splint to treat strong flexion synergy

High
L3915
Reimbursement

SKIN SAFE
Anti-microbial
material for skin
breakdown
prevention



COMFORT
Super Soft and
removable
Headliner cover is
machine
washable

MOVABLE
Adjustable
hinge at wrist
helps increase
wrist extension



FINGER CONE
Worn laterally with
fingers positioned
around cone. Adjusts
to contractures



Shown with optional finger separator

ADJUSTABLE
Patented spine
bends along
with movable
wrist joint

SIDE WINGS
Adjustable wings
prevent ulnar and
radial deviation

Ambidextrous

Size	Wrist to Tip to Middle Finger	Total Length
Adult	7" – 8"	12 ½"
Pediatric Large	4 ½" – 6"	9 ½"
Pediatric Medium	3" – 4 ½"	8"
Pediatric Small	2 ½" – 3 ½"	7"

