

Comfy™ Goniometer & Spring-Loaded Goniometer Elbow and Knee Orthoses

Patented goniometers allow for movement and progressive knee extension



COMFORT

Super Soft and removable Headliner or Terry cloth cover is machine washable. No exposed metal parts!

PROGRESSIVE EXTENSION

Gradual extension like serial casting with goniometer. Can also be set in any desired static position

ADJUSTABLE CUFFS

Patented malleable frame adjusts to individual needs without tools or heat

GONIOMETER DIAL

Set patented dial to desired range of flexion & extension from 0° to 180° in 10° increments. Supplied Allen wrench fits in splint cover.

SKIN SAFE

Anti-microbial and moisture wicking material for skin breakdown prevention

GENTLE PUSH

Optional Spring-loaded option provides a spring that gently pushes knee/elbow to selected degree of extension

	Joint crease up/down	Total Length	Comfy Elbow	Comfy Knee
Adult	7"	14"	(S)GE-101	(S)GK-101
Adult Sm	6"	12"	(S)GE-101-AS	(S)GK-101-AS
Ped. LG	5 ¼"	10 ½"	P(S)G-101-L	P(S)GK-101-L
Ped. Md	4 ½"	9"	P(S)GE-101-M	P(S)GK-101-M
Ped. Sm	3 ¾"	7 ½"	P(S)GE-101-S	P(S)GK-101-S

Watch our informational You Tube video online!



ComfySplints
We Make The World Comfy