



## Toilet Safety Frame

The **ProBasics Toilet Safety Frame** provides a safe, stable aid for individuals who struggle with the up and down demands of toileting. The frame is designed to accommodate both standard and elevated toilet seats. It is easy to clean and the aluminum frame mounts securely onto the bowl with an adjustable bracket. The toilet safety frame can safely support patient weights up to 250 pounds.

### FEATURES

- Frame is made of anodized, non-corrosive one-inch aluminum tubing
- Arms are height and width adjustable
- Versatile design with tool-free assembly quickly secures to most toilets
- Contoured waterfall armrests provide additional comfort and support while providing an excellent gripping surface
- Powder-coated aluminum bracket easily attaches frame to toilet
- HCPCS Code: E0243\*

### SPECIFICATIONS

Adjustable Arm Width:	20" – 22"
Adjustable Arm Height:	26" – 29.5"
Weight Capacity:	300 lb
Carton Dimensions:	19.5" x 18" x 10"
Shipping Weight:	15 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Description	UPC	UOM
BSTF	ProBasics Toilet Safety Frame	815067071821	4/cs