



Padded Transfer Bench

The **ProBasics Padded Transfer Bench** helps make transfers in and out of the bathtub safer and more comfortable. Constructed with an aluminum anodized frame and padded back and bench, it features a reversible back for right or left handed transfers. The angled legs have one-inch, easy-snap button adjustments and wide-based, non-skid feet for added sturdiness and safety. Tool-free assembly and limited lifetime warranty are standard.

Aluminum frame is lightweight, sturdy and corrosion-proof

- Vinyl seat and back are padded with 1.25" open-cell foam for added comfort
- Wide-based tips and suction cups add sturdiness and prevent slipping
- Legs have one-inch easy snap-button adjustments
- Converts easily for left or right hand transfers
- Quick tool-free assembly
- HCPCS Code: E0247*

Item #	Description	UPC	UOM
BSTBP	ProBasics Padded Transfer Bench	815067071982	2/cs