



Bath Bench

The **ProBasics Bath Bench** is designed for individuals who require seating support and stability when bathing or showering. The bench features a comfortable contoured textured seat that can support patient weights up to 300 pounds. The angled legs have one-inch, easy-snap button adjustments and wide-based, non-skid feet for added sturdiness and safety. Tool-free assembly and limited lifetime warranty are standard.

Aluminum frame is lightweight, durable and corrosion-proof

- Blow-molded bench has a textured surface and is contoured with drainage holes to reduce slipping
- Quick tool-free assembly
- Wide-based crutch tips help increase stability
- Two-hole handles in seat provide security when sitting or standing
- Legs are height adjustable with one-inch snap button adjustments
- HCPCS Code: E0245*

SPECIFICATIONS

| | |
|-------------------------|--|
| Overall Footprint: | 16.5"W x 19"D |
| Seat Dimensions: | 19.5"W x 11.5"D |
| Adjustable Seat Height: | 14.25" – 20.5" |
| Weight Capacity: | 300 lb |
| Carton Dimensions: | 23.5" x 17" x 19" |
| Shipping Weight: | 19 lb |
| Warranty: | Limited Lifetime on frame One year on all parts |

| Item # | Description | UPC | UOM |
|--------|----------------------|--------------|------|
| BSC | ProBasics Bath Bench | 815067071920 | 4/cs |