




ProBasics K1 Wheelchair

STANDARD WHEELCHAIR

The **ProBasics K1 Manual Wheelchair** combines style and durability at an affordable price. Standard features include flip back removable padded armrests, easy-to-clean vinyl upholstery, chart pocket and heel loops on the front riggings. Supporting patient weights up to 250 pounds, the ProBasics K1 Standard Wheelchair is a popular choice among providers. Available with swing-away footrests or elevating legrests.



Flip-back armrests



Easy-to-operate
wheel locks



Front riggings
with heel loop

FEATURES

- Slide tube steel frame in black powder coat finish
- Single axle positioning
- Easy-to-clean, comfortable vinyl upholstery resists mildew and bacteria
- Flip-back, padded, upholstered armrests are removable
- Chart pocket on back for additional convenience
- Includes tool-free swing-away footrests with heel loops or elevating padded legrests with composite footplates
- Elevating legrests are equipped with a pneumatic piston for infinite elevating positions
- Easy push-to-lock wheel locks
- Composite handrims with 24" solid urethane tires
- 8" x 1" solid rubber front casters
- Precision sealed wheel bearings in front and rear ensure long lasting performance and reliability
- HCPCS* Code: K0001

K0195 OR E0990 (elevating legrests only)

SPECIFICATIONS

Seat Width:	16" 18" 20"
Seat Depth:	16"
Seat-to-Floor Height ¹ :	19.5"
Back Height:	18"
Back Style:	Fixed
Overall Width w/Arms:	24 ½" 26 ½" 28 ½"
Overall Depth:	31 ¼" (without riggings)
Weight Capacity:	250 lb
Product Weight ² :	41 – 43 lb
Carton Dimensions:	34" x 12" x 41"
Shipping Weight:	45 – 50 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Seat	Arm Style	Leg Style	UPC	UOM
WC11616DS	16" x 16"	Flip Back Desk-Length	Swing-Away Footrests	815067070558	each
WC11816DS	18" x 16"	Flip Back Desk-Length	Swing-Away Footrests	815067070565	each
WC12016DS	20" x 16"	Flip Back Desk-Length	Swing-Away Footrests	815067070572	each
WC11616DE	16" x 16"	Flip Back Desk-Length	Elevating Legrests	815067070589	each
WC11816DE	18" x 16"	Flip Back Desk-Length	Elevating Legrests	815067070596	each
WC12016DE	20" x 16"	Flip Back Desk-Length	Elevating Legrests	815067070602	each

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.

¹NOTE: These heights can vary ¼-inch due to tire wear and manufacturing tolerances.

²NOTE: Weights based on 18" x 16" wheelchair with desk-length arms and no front riggings. Weights will vary depending on how the wheelchair is equipped.