



Folding Canes

ProBasics Folding Canes fold easily for storage or travel. The comfortable T-shape grip reduces hand stress and fatigue on the wrist for maximum comfort. Available in two attractive finishes with a limited lifetime warranty. Supports patient weights up to 300 pounds.



Folds compact, perfect for a purse or carry-on luggage



Contoured T-shape handle

FEATURES

- Folds in four convenient parts for easy storage
- Easy-to-use, one button height adjustment with a locking anti-rattle ring allowing for quick safe fitting
- Comfortable, contoured T-shape handle
- Easy-to-clean aluminum frame
- Adjusts in 1-inch increments
- Attractive powder-coated finish available in black or bronze
- HCPCS* Code: E0100

SPECIFICATIONS

Adjustment Range:	32.5" – 36.5"
Cane Diameter:	7/8"
Cane Tip:	3/4"
Product Weight:	13 oz
Weight Capacity:	300 lb
Carton Dimensions:	20.5" x 6.3" x 12.6"
Shipping Weight:	16.6 lb
Warranty:	Limited lifetime on frame 90 days on all parts

Item #	Description	Color	UPC	UOM
CNFBK	ProBasics Folding Cane	Black	815067071241	6cs
CNFBZ	ProBasics Folding Cane	Bronze	815067071258	6/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.

