

Fixed - Tilt

METRO



- **Contoured Upholstery**
Provides a closer fit for greater support and comfort
- **Adjustable Trunk Support**
Promotes superior lumbar positioning
- **Self-Tensio™**
Built in self-tensioning seating system increases sitting tolerance



Transit



Adheres to WC19 safety standards for safe transportation
(optional)



Mesh Upholstery Center



Increases ventilation to keep the occupant cooler

"Quick Release" Rear Wheels



Easy to remove for storage or transport

Two-Piece Push Handle



Quick and easy to fold, store, and transport

Lightweight • Compact Folding • Convenient

 **Convaid**
Now You're Going Places



Accessories



3-Pt Positioning Belt with Depth Adjustable Crotch Strap (Transit Required)



Travel Bag



Under Seat Storage Basket

Tire Options

- 12.5" Rear Solid Knobby Tire
- 7.5" x 2" Front 12.5" Rear Pneumatic Knobby Tire
- 8" x 2" Front Ribbed, 12.5" Rear Solid or Pneumatic Knobby Tire

Transit Options

- Q'Straint Transit Lap Belt
- Sure-Lok Transit Lap Belt

Measurement	Chair Model	Seat Width*	Overall Width	Overall Height	Overall Length	Seat Depth*	Seat Back Height	Seat to Footrest*	Seat to floor	Seat Angle	Headrest Extension	Weight Capacity Regular/Transit	Weight of Chair Δ
in.	MT 12	12	23.5	36.8	36.5	11	21	8 - 15	24	30°	6	75 / 66 lbs	23 lbs
cm	MT 12	30.5	60	93.3	92.7	28	53.3	20.3 - 38	61		15	34 / 30 kg	10.4 kg
in.	MT 14	14	25.5	37.5	34.5	13	23	8 - 14.5	24	30°	8	100 / 100 lbs	24 lbs
cm	MT 14	35.5	64.8	95.3	87.6	33	58.4	20.3 - 36.8	61		20	45.5 / 45.5kg	11 kg
in.	MT 16	16	27.5	38	35	14	25	13 - 18	23	30°	8	170 / 170 lbs	27 lbs
cm	MT 16	40.5	70	96.5	89	35.6	63.5	33 - 45.7	58.4		20	77 / 77 kg	12.3 kg
in.	MT 18	18	28.5	38	42	16	25	3 - 24	23.5	30°	8	250 / 170 lbs	29 lbs
cm	MT 18	45.5	72.4	96.5	106.7	40.6	63.5	7.3 - 61	59.7		20	114 / 77 kg	13.2 kg

in. - Inch, cm - Centimeter * Plus or Minus .5 inches (1.27 cm) due to fabric. Δ- Weight of chair is without front rigging

Upholstery Colors



Cherry Red



Royal Blue



Metro Transit