

TILT-IN-SPACE RODEO



- **Up to 45° Adjustable Tilt**
Provides pressure relief and upper body support
- **Customizable**
Custom-built to your specifications to meet each individual's physical needs
- **Multiple Upholstery Styles & Colors**
Removable, washable Cordura available in various colors
- **5 Years of Growth Capabilities**
Adjustable seat depth extends chair life for a growing child
- **Dynamic Positioning**
Accommodates functional or rest positions and increased sitting tolerance for maximum comfort



Positioning options & accessories sold separately

Transit




Meets or exceeds WC19 safety standards for safe transportation

90° to 110° Recline



Opens the seat-to-back angle to provide postural relief

Height Adjustable Flip-Up Armrest



Provide upper extremity support and positioning

Firm Laterals



Firm lateral thoracic support maintains midline positioning of the child's trunk

Lightweight • Compact Folding • Convenient

 **Convaid**
Now You're Going Places



Accessories

- Extended Headrest Cover with windows **A** (Canopy)
- Extended Headrest Cover without windows (Canopy)
- Headrest Cover (Canopy)
- Mosquito Net **B** (Canopy Required)
- Rain Cover (Canopy Required)
- Attendant Hand Brake **C**
- Upper Extremity Support Surface Tray **D** (Clear Tray) (Armrests Required)
- Utility Bag
- Oxygen Tank Bag
- Medical Necessity Under Seat Storage Basket **E**
- Under Seat Storage Basket
- Rear Anti-Tip Tubes **F**
- Caster Locks **G** (For 2" Wide Front Wheels)



Head Supports

- Headrest Extension (Required for Transit)
- Occi Headwings
- Padded Headwings

Trunk Positioning

- H-Harness with Padded Covers (Required for Transit)
 - Adjustable Firm Lateral Trunk Support Contoured (Swing-Away)
 - Firm Lateral Trunk Support Securement Strap (Scoliosis Strap)
 - Reducer Seat Insert (Shown)
 - Full Torso Support Vest
- (Lateral Shown with Scoliosis Strap)

Pelvic Positioning

- Support Cushion
- Align Cushion
- Position Cushion
- Solid Seat Insert (Required for all Cushions)
- Incontinence Cushion Liner (Solid Seat Insert & Cushion Required)

Frame Options

- Recline System (90° to 110°)
- Height Adjustable Flip-up Armrests
- Side Panels -Solid -Mesh

Transit Options

- Upper Extremity Support Surface (Foam Tray) (Armrests Required)
- Q'Straint Transit Lap Belt
- Sure-Lok Transit Lap Belt

Adjustable Tilt & Recline



Upholstery Colors

Cordura

- Panther Black
- Cherry Red
- Sassy Purple
- Princess Pink
- Royal Blue
- Forest Green
- Chocolate Brown
- Aqua Blue
- Tangerine Orange
- Apple Green

Foot & Leg Positioning

- Medial Thigh Support Abductor
- Lateral Thigh Support Adductor
- Calf Panel
- Foot Positioners (Required for Transit)
- 2" Legrest Extension
- Elevating Legrest

Compact Fold



Tire Options

- Front Tire**
- 7.5" x 2" Front Pneumatic ----- RD 10/12/14/16
- Rear Tire**
- 12.5" x 2.5" Rear Solid Knobby ----- RD 10/12/14/16
 - 12.5" x 2.5" Rear Pneumatic Knobby ----- RD 10/12/14/16