

Tilt-in-Space SAFARI



- **5° to 45° of Adjustable Tilt**
Provides pressure relief and upper body support
- **Customize**
Custom-built to your specifications to meet each individual's physical needs
- **Upholstery Colors**
Removable, washable Cordura available in various colors
- **5 Years of Growth Capabilities**
Adjustable seat depth extends chair life for a growing child
- **Adjustable Recline**
Opens the seat-to-back angle from 90° up to 135° to provide postural relief

Options:

Transit



Adheres to WC19 safety standards for safe transportation

Elevating Legrests



Allows for changes in position and promotes better circulation

Attendant Hand Brake



For increased safety and control

Adjustable Tilt



Provides pressure relief and upper body support

Support Cushion



Promotes stability and posture control for basic comfort and positioning

Standard Features:

Lightweight • Compact Folding • Convenient

Convaid
Now You're Going Places



Accessories



Headrest Cover
(Canopy)



Upper Extremity
Support Surface Tray
(Hardware Required)



Travel Bag
(For Wheelchair Only)



Foam Tray



Attendant
Hand Brakes



Utility Bag



Curb Tipper

Pelvis Positioning



3-Pt Positioning Belt
with Depth Adjustable
Crotch Strap
(Transit Required)



Seat Cushions

- Position
- Align

- Solid Seat Insert
- Incontinence Cushion Liner
(Cushion and Solid Seat Required)

Upholstery Colors

Cordura



Panther
Black



Cherry
Red



Sassy
Purple



Princess
Pink



Royal
Blue



Forest
Green



Chocolate
Brown



Aqua
Blue



Tangerine
Orange



Apple
Green

Foot & Leg Positioning



Medial Thigh Support
(Abductor)



Lateral Thigh Support
(Adductor)



Elevating Legrests
w/ SA Footplate
(Transit Required)



Calf Panel
• Foot Positioners
(Transit Required)

- Angle Adjustable Footplates
- Footplate Securement Strap

- Legrest Extension 3"

Head Supports

- Occi Headwings
- Padded Headwings
- Headrest Extension (Transit Required)

Tire Options

- 12.5" Rear Solid Knobby Tire
- 7.5" x 2" Front 12.5" Rear Pneumatic Knobby Tire
- 6" x 2" Front Solid Tire
- 6" x 2" Front Ribbed, 12.5" Rear Solid Knobby Tire

Frame Options

- Adjustable Recline 90° - 135°
- Anatomic Back
- Height Adjustable Flip-up Armrest
- Solid Side Panels

Trunk Positioning



H-Harness with
Padded Covers
(Transit Required)



Soft Adj. Lateral Support
Single or Double Flap



Full Torso Support Vest



Reducer Seat Insert

Transit Options

- Q'Straint Transit Lap Belt
- Sure-Lok Transit Lap Belt
- Upper Extremity Support Surface (Foam Tray)

Measurement	Chair Model	Dimensions (inches)											Capacity		Weight of Chair Δ		
		Seat Width*	Overall Width	Overall Height	Overall Length	Seat Depth*	Seat Back Height	Seat To Footplate*	Seat to Floor	Armrest to Seat	Seat to Back Angle	Recline Adjustment	Tilt Adjustment (from vertical)	Headrest Extension		Shoulder Strap Height	Weight Capacity Regular/Transit
in.	SF 12	12	21.5	44.5	36.5	9 - 14	23	11 - 13	23.5	8 - 11	90° - 95°	3°, 22°, 45°	5° - 45°	8	8.5 - 20	75 / 66 lbs	34 lbs
cm	SF 12	30.5	54.6	113	92.7	23 - 35.6	58	27.9 - 33	59.7	20.3 - 27.9				20	22 - 51	34 / 30 kg	15.5 kg
in.	SF 14	14	24.3	45.5	39	11 - 16	22	11 - 13	26	8 - 10.5	90° - 95°	1°, 19°, 40°	5° - 45°	8	9.5 - 19.5	100 / 100 lbs	37 lbs
cm	SF 14	35.5	61.6	115.6	99	27.9 - 40.6	56	27.9 - 33	66	20.3 - 26				20	24 - 50	45.5 / 45.5 kg	16.8 kg
in.	SF 16	16	25.5	45.5	42.5	13 - 18	22	12.5 - 14.5	26	8 - 11	90° - 95°	1°, 19°, 40°	5° - 45°	8	10 - 20	170 / 170 lbs	41 lbs
cm	SF 16	40.5	64.8	115.6	108	33 - 45.7	56	31.8 - 36	66	20.3 - 27.9				20	25 - 51	77 / 77 kg	18.6 kg
in.	SF 18	18	26	44.5	41.5	15 - 20	21.5	10 - 16	23	7.5 - 10	90° - 95°		5° - 45°	8	12.5 - 18.5	250 / 170 lbs	45 lbs
cm	SF 18	45.5	66	113	105.5	38 - 51	55	25 - 41	58	19 - 25				20	32 - 47	114 / 77 kg	20.5 kg

in. - Inch, cm - Centimeter * Plus or Minus .5 inches (1.27 cm) due to fabric. Δ- Weight of chair is without front rigging