

Convoid Rodeo Tilt - HCPCS E1232/E1161

Retail Price List / Order Form

Account #: _____ Purchase Order #: _____ Date: _____
 Contact: _____ Mark For: _____

| | |
|-----------------|-----------------|
| Bill To: | Ship To: |
| Address: | Address: |
| City/State/Zip: | |
| Phone: | City/State/Zip: |
| Email: | Phone: |

Rodeo Tilt Standard Features:

- 5° to 45° of Adjustable Tilt
- Silver Powder Coated Frame
- One Piece Height Adjustable Push Handle
- One Piece Angle Adjustable Footplate
- 7.5 x 2" Front, 11.5 x 2.5" Rear Solid Tire
- Quick Release Rolko Rear Tires
- Single Foot Activated Wheel Lock Bar
- Zippered Pockets in Seat Bottom Upholstery
- 2-Point Positioning Belt
- Seat Depth Growth of 5"
- Elevating Legrest
- HCPCS code E1161 is only valid on size 16 (RD16)
- Tool-less Upholstery Removal
- Made in the USA

Standard Models:

| | <u>Price</u> |
|--|----------------|
| <input type="checkbox"/> Rodeo RD10 - 904340 | \$2,468 |
| <input type="checkbox"/> Rodeo RD12 - 902689 | \$2,468 |
| <input type="checkbox"/> Rodeo RD14 - 902691 | \$2,468 |
| <input type="checkbox"/> Rodeo RD16 - 903134 | \$2,468 |

Upholstery Colors

| | <u>Price</u> |
|--|--------------|
| <input type="checkbox"/> Royal Blue | \$0 |
| <input type="checkbox"/> Forest Green | \$0 |
| <input type="checkbox"/> Cherry Red | \$0 |
| <input type="checkbox"/> Sassy Purple | \$0 |
| <input type="checkbox"/> Princess Pink | \$0 |
| <input type="checkbox"/> Panther Black | \$0 |
| <input type="checkbox"/> Chocolate Brown | \$0 |
| <input type="checkbox"/> Aqua Blue | \$0 |
| <input type="checkbox"/> Apple Green | \$0 |

| RD10 | RD12 | RD14 | RD16 |
|--------|--------|--------|--------|
| 903850 | 903850 | 903850 | 903850 |
| 903857 | 903857 | 903857 | 903857 |
| 903855 | 903855 | 903855 | 903855 |
| 903851 | 903851 | 903851 | 903851 |
| 903759 | 903759 | 903759 | 903759 |
| 903854 | 903854 | 903854 | 903854 |
| 904189 | 904189 | 904189 | 904189 |
| 904753 | 904753 | 904753 | 904753 |
| 904754 | 904754 | 904754 | 904754 |

Transit Required Items

| | <u>Price</u> | <u>HCPCS</u> |
|--|--------------|--------------|
| <input type="checkbox"/> Transit Option (Manufacturer Installation Required) | \$200 | |
| <input type="checkbox"/> Headrest Extension | \$130 | E0966 |
| <input type="checkbox"/> H-Harness with Padded Covers | \$143 | E0960 |
| <input type="checkbox"/> 3-Pt Positioning Belt | \$49 | E0978 |
| Foot Positioners - Select Size Below: | | |
| <input type="checkbox"/> Foot Positioners - Small 7" - 10" (pair) | \$156 | K0108 |
| <input type="checkbox"/> Foot Positioners - Medium 8" - 11" (pair) | \$156 | K0108 |
| <input type="checkbox"/> Foot Positioners - Large 9" - 12" (pair) | \$156 | K0108 |

| RD10 | RD12 | RD14 | RD16 |
|--------|--------|--------|--------|
| 904341 | 903167 | 903168 | 903169 |
| 903332 | 902863 | 900484 | 900138 |
| 903988 | 903988 | 903988 | 903988 |
| 904017 | 904017 | 904017 | 904017 |
| | | | |
| 903415 | 903415 | 903415 | 903415 |
| 903416 | 903416 | 903416 | 903416 |
| 903417 | 903417 | 903417 | 903417 |

Transit Options

| | <u>Price</u> | <u>HCPCS</u> |
|---|--------------|--------------|
| <input type="checkbox"/> Q'Straint Transit Lap Belt | \$90 | E0978 |
| <input type="checkbox"/> Sure-Lok Transit Lap Belt | \$90 | E0978 |
| <input type="checkbox"/> Upper Extremity Support Surface (Foam Tray) Small ^ | \$151 | E0950 |
| <input type="checkbox"/> Upper Extremity Support Surface (Foam Tray) Medium ^ | \$151 | E0950 |

^ Crash tested Foam Tray with vinyl cover. Flip-Up Armrests required.

| RD10 | RD12 | RD14 | RD16 |
|--------|--------|--------|--------|
| 903094 | 903094 | 903094 | 903094 |
| 903232 | 903232 | 903232 | 903232 |
| 903669 | 903669 | | |
| | | 903670 | 903670 |

Pelvic Positioning

| | <u>Price</u> | <u>HCPCS</u> |
|--|--------------|--------------|
| <input type="checkbox"/> 3-Pt Positioning Belt | \$49 | E0978 |
| <input type="checkbox"/> Support - General Use Cushion * ¹ | \$241 | E2601 |
| <input type="checkbox"/> Position Cushion - Medial Thigh Support/Anti Thrust * ¹ | \$466 | E2607 |
| <input type="checkbox"/> Align Cushion - Lateral Pelvis/Lateral Thigh Support * ¹ | \$466 | E2607 |
| <input type="checkbox"/> Incontinence Cushion Liner (Cushion req) ¹ | \$72 | E2619 |

* Solid Seat Insert included with all Cushions.
¹ Manufacturer Installation Required.

| RD10 | RD12 | RD14 | RD16 |
|--------|--------|--------|--------|
| 904017 | 904017 | 904017 | 904017 |
| 904347 | 904350 | 904353 | 904356 |
| 904348 | 904351 | 904354 | 904357 |
| 904349 | 904352 | 904355 | 904358 |
| 902324 | 902324 | 902324 | 902324 |

Foot & Leg Positioning

| | <u>Price</u> | <u>HCPCS</u> |
|--|--------------|--------------|
| <input type="checkbox"/> Individual Angle Adjustable Footplates ¹ | \$179 | K0040 |
| <input type="checkbox"/> Full Padded Footbox (One-Piece Footplate Required) | \$325 | K0108 |
| <input type="checkbox"/> Foot Positioners - Small 7" - 10" (pair) | \$156 | K0108 |
| <input type="checkbox"/> Foot Positioners - Medium 8" - 11" (pair) | \$156 | K0108 |
| <input type="checkbox"/> Foot Positioners - Large 9" - 12" (pair) | \$156 | K0108 |
| <input type="checkbox"/> Legrest Extension - 2" * ¹ | \$55 | |
| <input type="checkbox"/> Medial Thigh Support (Abductor) | \$162 | E0957 |
| <input type="checkbox"/> Lateral Thigh Support (Adductor) | \$162 | E0956 |
| <input type="checkbox"/> Calf Panel | \$61 | E0995 |
| <input type="checkbox"/> Padded Elevating Legrest Cover | \$10 | |

| RD10 | RD12 | RD14 | RD16 |
|--------|--------|--------|--------|
| | 904957 | 904958 | 904959 |
| 904953 | 904954 | 904955 | 904956 |
| 903415 | 903415 | 903415 | 903415 |
| 903416 | 903416 | 903416 | 903416 |
| 903417 | 903417 | 903417 | 903417 |
| 904095 | 904095 | 904096 | 904097 |
| 904359 | 902667 | 902647 | 903140 |
| 904360 | 902668 | 902648 | 903141 |
| 904361 | 904361 | 904361 | 904361 |
| 904619 | 904619 | 904619 | 904619 |

* 2" Legrest Extension interferes with caster swivel when seat is in upright position.
 Seat must be tilted back at least 5 degrees OR Legrests must be elevated to allow caster swivel.
¹ Manufacturer Installation Required.

| <u>Trunk Positioning</u> | <u>Price</u> | <u>HCPCS</u> |
|--|--------------|--------------|
| <input type="checkbox"/> H-Harness with Padded Covers (3-Pt Positioning Belt Req) | \$143 | E0960 |
| Right Adj Firm Lateral Trunk Sup - Contoured (Hrdw Req) | | |
| <input type="checkbox"/> X Small= 2.75"x 3.75" * | \$118 | E0956 |
| <input type="checkbox"/> Small = 3"x 5" | \$118 | E0956 |
| <input type="checkbox"/> Medium = 3.75"x 5.5" | \$118 | E0956 |
| <input type="checkbox"/> Large = 4.75"x 6" | \$118 | E0956 |
| Left Adj Firm Lateral Trunk Sup - Contoured (Hrdw Req) | | |
| <input type="checkbox"/> X Small= 2.75"x 3.75" * | \$118 | E0956 |
| <input type="checkbox"/> Small = 3"x 5" | \$118 | E0956 |
| <input type="checkbox"/> Medium = 3.75"x 5.5" | \$118 | E0956 |
| <input type="checkbox"/> Large = 4.75"x 6" | \$118 | E0956 |
| <input type="checkbox"/> Right S/A Adj. Hardware for Lateral Trunk Support | \$238 | E1028 |
| <input type="checkbox"/> Left S/A Adj. Hardware for Lateral Trunk Support | \$238 | E1028 |
| <input type="checkbox"/> Full Torso Support Vest | \$238 | |
| <input type="checkbox"/> Reducer Seat Insert (Reduces Seat Width 2" & Seat Depth 1") | \$353 | |

* RD10 requires X Small laterals.

| <u>Head Supports</u> | <u>Price</u> | <u>HCPCS</u> |
|---|--------------|--------------|
| <input type="checkbox"/> Padded Headwings | \$244 | E0955 |
| <input type="checkbox"/> Occi Headwings | \$228 | E0955 |
| <input type="checkbox"/> Headrest Extension | \$130 | E0966 |

| <u>Frame Option</u> | <u>Price</u> | <u>HCPCS</u> |
|--|--------------|--------------|
| <input type="checkbox"/> Recline Option (Requires Manufacturer Installation) | \$461 | E1225 |
| <input type="checkbox"/> Height Adjustable Flip-Up Armrests (pair) | \$228 | E0973 |
| <input type="checkbox"/> Solid Side Panels | \$211 | |
| <input type="checkbox"/> Mesh Side Panels | \$148 | |

| <u>Tire Options</u> | <u>Price</u> | <u>HCPCS</u> |
|--|--------------|--------------|
| <input type="checkbox"/> 12.5" Rear Solid Knobby Tire | \$199 | |
| <input type="checkbox"/> 7.5" x 2" Front, 12.5" Rear Pneumatic Knobby Tire | \$299 | |

| <u>Accessories</u> | <u>Price</u> | <u>HCPCS</u> |
|--|--------------|--------------|
| <input type="checkbox"/> Headrest Cover (Canopy) | \$228 | |
| <input type="checkbox"/> Extended Headrest Cover - No Windows (Canopy) | \$250 | |
| <input type="checkbox"/> Extended Headrest Cover - with Windows (Canopy) | \$295 | |
| <input type="checkbox"/> Upper Extremity Support Surface (Clear Tray) ^ | \$111 | E0950 |
| <input type="checkbox"/> Adj. Removable Hardware (Required for Clear Tray) | \$227 | E1028 |
| <input type="checkbox"/> Medical Necessity Storage Basket | \$295 | |
| <input type="checkbox"/> Under Seat Storage Basket - Tilt-In-Space | \$75 | |
| <input type="checkbox"/> IV Pole (height adjustable and collapsible) | \$150 | K0105 |
| <input type="checkbox"/> Oxygen Tank Bag | \$179 | E2208 |
| <input type="checkbox"/> Utility Bag | \$38 | |
| <input type="checkbox"/> Attendant Hand Brakes (pair) ¹ | \$236 | |
| <input type="checkbox"/> Curb Tipper | \$98 | |
| <input type="checkbox"/> Caster Lock for 2" Wide Front Tires ¹ | \$94 | K0073 |
| <input type="checkbox"/> Rear Anti-Tip Tubes (Pair)* | \$115 | E0971 |
| <input type="checkbox"/> Rain Cover Size 1 (Canopy Req) | \$150 | |
| <input type="checkbox"/> Rain Cover Size 2 (Canopy Req) (RD14 Extended Canopy Req) | \$150 | |
| <input type="checkbox"/> Rain Cover Size 3 (Extended Canopy Req) | \$150 | |
| <input type="checkbox"/> Mosquito Net Size 1 (Canopy req) | \$150 | |
| <input type="checkbox"/> Mosquito Net Size 2 (Canopy req) (RD14 Extended Canopy Req) | \$150 | |
| <input type="checkbox"/> Mosquito Net Size 3 (Extended Canopy req) | \$150 | |

^ Flip-Up Armrests and Hardware required.

*Highly recommended for users weighing 50lbs or less.

¹ Manufacturer Installation Required.

| RD10 | RD12 | RD14 | RD16 |
|--------|--------|--------|--------|
| 903988 | 903988 | 903988 | 903988 |
| 904385 | 904385 | 904385 | 904385 |
| | 904071 | 904071 | 904071 |
| | 904072 | 904072 | 904072 |
| | 904073 | 904073 | 904073 |
| 904465 | 904465 | 904465 | 904465 |
| | 904466 | 904466 | 904466 |
| | 904467 | 904467 | 904467 |
| | 904468 | 904468 | 904468 |
| 904362 | 904362 | 904362 | 904362 |
| 904363 | 904363 | 904363 | 904363 |
| 904364 | 903261 | 903260 | 903259 |
| 904211 | 904211 | 904110 | 904110 |

| RD10 | RD12 | RD14 | RD16 |
|--------|--------|--------|--------|
| 904365 | 902673 | 902653 | 903147 |
| 905270 | 905270 | 905271 | 905271 |
| 903332 | 902863 | 900484 | 900138 |

| RD10 | RD12 | RD14 | RD16 |
|--------|--------|--------|--------|
| 904342 | 904343 | 904344 | 904345 |
| 904346 | 904346 | 904346 | 904346 |
| 904388 | 902672 | 902652 | 903145 |
| 904389 | 902680 | 902681 | 903146 |

| RD10 | RD12 | RD14 | RD16 |
|--------|--------|--------|--------|
| 904117 | 904117 | 904117 | 904117 |
| 904118 | 904118 | 904118 | 904118 |

| RD10 | RD12 | RD14 | RD16 |
|--------|--------|--------|--------|
| 904366 | 902679 | 902659 | 903152 |
| 904728 | 904729 | 904730 | 904731 |
| 904744 | 904745 | 904746 | 904747 |
| 904367 | 903405 | 903401 | 903402 |
| 903369 | 903369 | 903369 | 903369 |
| 904476 | 904476 | 904477 | 904477 |
| 904529 | 904529 | 904530 | 904530 |
| 904566 | 904566 | 904566 | 904566 |
| 902870 | 902870 | 902870 | 902870 |
| 904531 | 904531 | 904531 | 904531 |
| | | 904370 | 904371 |
| 905310 | 905310 | 905310 | 905310 |
| 904650 | 904650 | 904650 | 904650 |
| 904649 | 904649 | 904649 | 904649 |
| 904939 | 904939 | 904939 | |
| | | 904940 | 904940 |
| | | | 905283 |
| 904655 | 904655 | 904655 | |
| | | 904656 | 904656 |
| | | | 905319 |

The Rodeo Transit option includes four manufacturer installed red wheelchair transportation anchors. The Rodeo Transit model has been crash tested and conforms to RESNA WC19 standards. Anthropomorphic Test Device weights are RD10T, RD12T 66lbs/30kg, RD14T 100lbs/45kg, and RD16T 140lbs/64kg.

The Convaid Rodeo mobility base was tested with Rodeo seating module equipped with head support, H-Harness with padded Covers, 3-point restraint system, foot positioners, and 4 transit anchor brackets.

It is Convaid's position that the wheelchair user should be transferred from the wheelchair into appropriate vehicle seating whenever possible. Please note that at this time, the Department of Transportation has not officially approved any tie-down system for a user while in a wheelchair, therefore Convaid cannot recommend a transit system.

The WC19 standard is specially designed for the wheelchair to be secured in a forward facing configuration and not in any other direction. For more information on RESNA WC-19 and WC-20 standards please visit <http://www.nercwtts.pitt.edu/WC19.html>.

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

Standard Domestic Shipping Rate: \$40

Prices in USD & subject to change
Effective: 4/20/2020

OF0100-ENG-REV.01