

CORFIT SYSTEM® INDUSTRIAL LS BACK SUPPORT LUMBOSACRAL SPINAL SUPPORT

- Ideal for lifting and industrial activities
- Internal suspender holds support in position
- Proportional sizing system for a proper fit
- Easy closure assists provide leverage with minimal effort
- Helps stabilize the spine and promote proper posture
- Flexible stays help increase support and stability
- Adjustable elastic pulls for extra compression
- Contours to the waist and hips for comfort



Available From:



STABILIZE THE SPINE AND EASE LOWER BACK PAIN

Ideal for lifting and industrial activities, the CorFit System® Industrial LS Back Support provides compression to help stabilize the spine and ease lower back pain by creating intra-abdominal pressure. An internal suspender holds this support in position and does not present a hazard around machinery or equipment like over the shoulder suspenders. It also helps improve posture and serves as a reminder to use proper lifting techniques by limiting the motion of muscles and ligaments.

PATENTED, PROPORTIONAL SIZING SYSTEM

Our patented, proportional sizing system ensures a proper fit in height and length to accommodate various body sizes. When ready to lift or use lumbosacral muscles, simply place hands through closure assists, tighten, and secure.

COMFORTABLE SUPPORT AND COMPRESSION

Easy closure assists provide additional leverage with minimal effort for extra compression and a secure fit; ideal for users with arthritis or weakened muscles. Elastic side pulls with convenient thumb loops allow extra compression for desired support. Flexible internal stays offer support and stability as well as help prevent the sides from rolling. The soft elastic CorEdge® finish on top and bottom edges offers added comfort and flexibility.

SPECIFICATIONS

Product Selection: CorFit System®
Industrial LS Back Support
LSB-7500, Black

Product Weight:
Small: .35 lbs | Medium/Large: .5 lbs |
XLarge: .6 lbs | 2XLarge: .7 lbs |
3XLarge: .9 lbs | 4XLarge: 1 lb

Dimensions:
Small: 6" x 35" | Medium/Large: 8" x 41" |
XLarge: 9" x 46" | 2XLarge: 10" x 52" |
3XLarge: 12" x 58" | 4XLarge: 14" x 64"

Materials: Elastic, heavy loop laminate

Patent Number: 5,785,672

Case Quantity: 12

Case Dimensions: 11" x 14" x 13"

Case Weight:
Small: 7 lbs | Medium: 9 lbs |
Large/XLarge: 12 lbs | 2XLarge: 12 lbs |
3XLarge: 13 lbs | 4XLarge: 15 lbs

Latex: This product is not made with natural rubber latex

Country of Origin: Made in USA