

SWEDE-O® ELASTIC ELBOW SLEEVE

- Light support for mild injuries, aching and/or weakened elbows
- Seamless design for a comfortable fit
- Gentle compression helps reduce inflammation
- Soft, ventilated knit material is suitable for extended use
- Contours to the natural shape of the elbow for comfort
- Bilateral design fits right or left elbow



Available From:



COMFORTABLE ELBOW SUPPORT

The Swede-O® Elastic Elbow Sleeve provides light support for mild injuries, aching and/or weakened elbows. Gentle compression helps reduce inflammation. The seamless, ventilated material contours to the natural shape of the elbow for comfort during extended use. Discreet under clothing; bilateral design fits right or left elbow.

SPECIFICATIONS

Product Selection: Swede-O® Elastic Elbow Sleeve ELB-6517, Beige

Product Weight: .3 lbs

Dimensions:

Small: 9" x 4.5"

Medium: 9" x 4.875"

Large: 9" x 5.25"

XLarge: 9" x 5.5"

Materials: Polyester, Rubber, Nylon

Case Quantity: 12

Case Dimensions: 12" x 8" x 10"

Case Weight: 5 lbs

Latex: Caution: This product contains natural rubber latex which may cause allergic reactions

Country of Origin: Made in Taiwan