

SWEDE-O® TRIM LOK® ANKLE BRACE

- Breathable mesh fabric helps keep foot and ankle cool, dry, and comfortable
- Lightweight and low profile
- Stainless steel spiral stays provide support
- Full elastic back reduces irritation and potential for blistering
- Bilateral design fits right or left ankle
- Industry leading 1 year warranty, satisfaction guarantee



Available From:



LIGHTWEIGHT, COOL & COMFORTABLE

The Swede-O® Trim Lok® brace is 25% lighter than other standard ankle braces. This brace is made from nylon mesh material to provide airflow, which helps to keep the foot and ankle cool, dry, and comfortable.

The built-in stainless steel spiral stays provide support to help prevent the ankle from rolling. The lightweight, low profile design fits comfortable in most shoes. A full elastic back is comfortable against the Achilles tendon and helps reduce irritation and potential for blistering. Backed by an industry leading 1 year warranty.

INDICATIONS FOR USE

Use the Trim Lok after injury or to protect your ankle from future injury. This ankle brace helps prevent and relieve aches, pains, sprains, and other ankle conditions.

SPECIFICATIONS

Product Selection: Swede-O® Trim Lok® AKL-6333, Black

Product Weight: .35 lbs

Dimensions:

XSmall: 7.375" x 8"

Small: 8.75" x 8.75"

Medium: 9" x 9"

Large: 9" x 9.5"

XLarge: 9" x 9.5"

Materials: Vinyl, polyester, nylon, polyurethane foam, cotton, steel, brass

Case Quantity: 12

Case Dimensions: 11" x 14" x 13"

Case Weight: 6 lbs

Latex: This product is not made with natural rubber latex

Country of Origin: Made in USA