

## SWEDE-O® WALKING BOOTS

- Helps immobilize, support, and protect the ankle and foot
- Rocker bottom promotes natural gait
- Cushioned sole for shock absorption
- Wide footbed and deluxe liner provides a comfortable fit
- Lightweight yet durable construction
- Uprights are anatomically contoured for a better fit
- Bilateral design fits right or left foot



### STUDY REFERENCE

<sup>1</sup> <https://www.ncbi.nlm.nih.gov/pmc/articles/PMC5721918/>

Available From:



Short

Tall

### HELPS SUPPORT ANKLE & FOOT

Swede-O® Walking Boots help immobilize the ankle and is ideal for severe ankle sprains, stable ankle fractures, post surgery, and post cast recovery. Studies<sup>1</sup> show that bearing weight on a stable injury helps facilitate recovery, allowing the patient to return to normal activity at a faster rate.

Designed with a non slip rocker bottom and low profile heel to promote natural gait. Shock absorbing, closed cell foam insole and heel cushion absorb and dissipate the impact of heel strike when walking. A wide footbed provides increased stability and accommodates for edema and bandage while a deluxe foam liner provides comfort.

Swede-O Walking Boot is durable, lightweight, and includes fitting pads to accommodate extra space around the ankle, heel, or leg. Anatomically contoured uprights provide a comfortable fit.

### SPECIFICATIONS

**Product Selection: Swede-O® Walking Boot, Black**

Short: UTL-1130 | Tall: UTL-1131

**Product Weight:**

Short: 2.3 lbs | Tall: 2.8 lbs

**Short Dimensions:** Small: 10.5" x 11" | Medium: 11.5" x 11.5" | Large: 13" x 13"

**Tall Dimensions:** Small: 15.5" x 10.5" | Medium: 17.5" x 11.5" | Large: 18.5" x 13"

**Materials:** pvc vinyl, polyethylene

**Case Quantity:** 6

**Case Dimensions:**

Short: 17" x 14" x 18" | Tall: 23" x 14" x 18"

**Case Weight:**

Short: 15 lbs | Tall: 18 lbs

**Latex:** This product is not made with natural rubber latex

**Country of Origin:** Made in China